Basic principles of SPE


Solid phase extraction (SPE) is a powerful method for sample preparation and is used by most chromatographers today.

More than 20 years ago MACHEREY-NAGEL designed and introduced CHROMABOND® SPE cartridges containing silica-based adsorbents. Since then we developed the widest range of phases and products for SPE based on silica and polymeric materials.

SPE has capabilities in a broad range of applications:

- environmental analyses
- pharmaceutical and biochemical analyses
- organic chemistry
- food analysis

Solid phase extraction


SPE is a form of digital (step-wise) chromatography designed to extract, partition, and/or adsorb one or more components from a liquid phase (sample) onto a stationary phase (adsorbent or resin). An adsorbed substance can be removed from the adsorbent by step-wise increase of elution strength of the eluent (step gradient technique). SPE extends a chromatographic system's lifetime, improves qualitative and quantitative analysis, and the demand placed on an analytical instrument is considerably lessened.


In general, SPE is used for three important purposes in state-of-the-art analyses:

- concentration of the analyte (up to factor 10.000 - increase of chromatographic sensibility / improved limits of detection)
- removal of interfering compounds (protection of subsequent analyses like HPLC, GC, TLC, UV or IR spectroscopy, ...)
- changing an analyte's environment to a simpler matrix more suitable for subsequent analyses

Advantages of SPE compared to classical liquid-liquid extraction:

- lower consumption of solvents
- faster - enormous time savings
- lower costs per sample
- potential for automation
- high consistency in individual sample handling
- more specific selectivity because of the broad range of adsorbents and different retention mechanisms
- optimisation of extraction by variation or adjusting of the solid phase and chromatographic conditions


Basic principles and summary of MN phases for SPE

Design of columns and cartridges

All CHROMABOND® columns and CHROMADIX® cartridges are manufactured from polypropylene (PP) with lowest content of extractables (plasticizers, stabilisers, ...) offering blank value free results by usage of most common solvents. The high quality CHROMABOND® adsorbents are kept in place by chemically very inert polyethylene filter elements (PE, standard pore size 20 µm). Funnel-shaped large volume columns, 96-well plates and cartridges for automated and on-line SPE packed with CHROMABOND® adsorbents are available on request.

CHROMABOND® polypropylene columns

- PP columns with PE filter elements
- different sizes from 1, 3, 6 up to 150 ml
- adsorbent weights from 20 mg to 50 g
- male luer tip as exit
- compatible with most robots (e.g., Gilson ASPEC™, Caliper AutoTrace®)

CHROMABOND® glass columns are available on request.

CHROMAFIX® cartridges

- PP cartridges with PE filter elements
- three different sizes with different adsorbent weights: **Small** (0.4 ml), **Medium** (0.8 ml), **Large** (1.8 ml)
- female Luer tip at the inlet, male Luer tip as exit
- offers alternative way of handling using positive pressure by syringes or peristaltic pumps
- especially suited for convenient SPE of small sample volumes

Solid phase extraction

Phase	Matrix	Modification / Application	Similar phases*
RP phases			
HR-X	PS/DVB		ENVI-Chrom P, Strata™-X, Oasis® HLB, Nexus
PS-RP	PS/DVB	removal of organic components	Strata™ SDB-L, Bond Elut® ENV, Bond Elut® LMS, DCS-PS/DVB, ENV PS-DVB, Bakerbond™ H ₂ O-phobic DVB, Isolute® 101
C ₁₈	silica	octadecyl, not endcapped	Strata™ C18-U, Accubond® C18, Bakerbond™ PolarPlus, Isolute® C18, LiChrolut® RP-18
C ₁₈ Hydra	silica	octadecyl, not endcapped, for polar analytes	
Normal phases			
SiOH	silica	unmodified	Strata™ Si-1, Bond Elut® silica, DSC-Si, LC-Si, CLEAN-UP® silica, Accubond® silica, Bakerbond™ silica gel, Isolute® silica, LiChrolut® Si
NH ₂	silica	aminopropyl	Strata™ NH ₂ , Sep-Pak® NH ₂ , Bond Elut NH ₂ , DSC-NH ₂ , LC-NH ₂ , CLEAN-UP® aminopropyl, Accubond® NH ₂ , Bakerbond™ amino, Isolute® NH ₂ , LiChrolut® NH ₂
OH	silica	diol	DSC-Diol, LC-Diol, Accubond® Diol (OH)
Alox A/N/B	aluminium oxide	acidic / neutral / basic	LC-Alumina-A/N/B, Accubond® aluminium oxide A/N/B
Florisil®	magnesium silicate		Strata™ FL-PR, Sep-Pak® Florisil®, Bond Elut® Florisil®, ENVI-Florisil®, LC-Florisil®, CLEAN-UP® Florisil®, Accubond® Florisil®, Bakerbond™ Florisil®, Isolute® FL, LiChrolut® Florisil®
PA	polyamide 6		DPA-6S
Ion exchangers			
SA	silica	benzenesulphonic acid cation exchanger (SCX)	Strata™ SCX, Bond Elut® SCX, DSC-SCX, LC-SCX, CLEAN-UP® Benzenesulfonic Acid, Accubond® SCX, Bakerbond™ Aromatic Sulfonic Acid, Isolute® SCX, LiChrolut® SCX
PSA	silica	propylsulphonic acid cation exchanger	
PS-OH ⁻	PS/DVB	strong anion exchanger, OH ⁻ form	Oasis® MAX
PS-H ⁺	PS/DVB	strong cation exchanger, H ⁺ form	Oasis® MCX, Strata™ X-C
PS-Ag ⁺	PS/DVB	strong cation exchanger, Ag ⁺ form	
PS-Ba ²⁺	PS/DVB	strong cation exchanger, Ba ²⁺ form	
Phases for special applications			
Drug	silica	bifunctional C ₈ /SA, for enrichment of drugs from urine	Strata™ Screen-C, Bond Elut® Certify I, DSC-MCAX, Clean Screen® DAU, Accubond® Evidex, Bakerbond™ Narc-2, Isolute® HCX, LiChrolut® TSC
Drug II	silica	bifunctional C ₈ /SB, for extraction of THC and derivatives as well as acidic analytes from biological fluids	Strata™ Screen-A, Bond Elut Certify II, Clean Screen® THC, Bakerbond™ Narc-1, Isolute® HAX
CN/SiOH	silica	combination phase for PAHs from soil	
NH ₂ /C ₁₈	silica	combination phase for PAHs from water	
Na ₂ SO ₄ /Florisil®	silica	combination phase for extraction of hydrocarbons from water (DIN H-53 / ISO DIS 9377-4)	
SA/SiOH	silica	combination phase for PCBs from waste oil	
Diamino	silica	primary and secondary amine functions (PSA), for de-termination of pesticides in food (QuEChERS method)	Supelclean PSA, Bond Elut PSA

* phases which provide a similar selectivity based on chemical or physical properties (list not complete)

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
Standard protocol for method development with CHROMABOND® HR-X

Column type: CHROMABOND® HR-X
3 ml / 200 mg


Sample pretreatment: if necessary, adjust pH value
Column conditioning: 5 ml methanol
Equilibration: 5 ml water
Sample application: slowly aspirate the sample through the column
Column washing: 5 ml water - methanol (95:5, v/v)
Elution: after drying 3 x 2 ml methanol

Further analysis: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC

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
2

PAHs from water containing humic acids

Column type:
CHROMABOND® NH₂/C₁₈, 6 ml, 500 mg/1 g
glass column

Sample pretreatment:

mix 500 ml water sample with 25 ml 2-propanol

Column conditioning: 10 ml dichloromethane, 10 ml methanol, then 10 ml dist. water – 2-propanol (9:1, v/v)

Sample application: aspirate 500 ml of the pretreated water sample through the column (~ 5 ml/min)

Washing: 2 ml dist. water – 2-propanol (9:1, v/v), then dry column (about 20 min, vacuum)

Elution: 4 x 0.5 ml CH₂Cl₂ (percolate first 0.5 ml into the column packing without vacuum, then apply light vacuum), if necessary evaporate in a stream of nitrogen and fill up with a suitable solvent

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
PAHs from soil

Column type:
CHROMABOND® CN/SiOH,
6 ml, 500/1000 mg

Sample pretreatment: dry 30 g soil with sodium sulphate and reflux 4 h with 250 ml petroleum ether in a Soxhlet extractor. For low PAH contents (colourless or weakly coloured extracts) concentrate extract to 1/10 of its volume in a rotation evaporator.

Column conditioning: 4 ml petroleum ether

MN Appl. Nr. 301310

Sample application:

aspirate 20 ml of the extract through the column


Washing: 2 ml petroleum ether

Elution: 2 x 2 ml acetonitrile / toluene (3:1, v/v), then evaporate or fill to the volume required

Further analysis: HPLC, e.g. with column 250 x 3 mm NUCLEOSIL® 5 C₁₈ PAH

recovery rates see application 301310 at www.mn-net.com

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Hydrocarbons from water

Column type:
CHROMABOND® Na₂SO₄/Florisil®,
2000/2000 mg, 6 ml glass column

Internal standard solution: dissolve 20 mg *n*-tetracontane (C₄₀H₈₂) in petroleum ether, add 20 ml *n*-decane (C₁₀H₂₂) and fill up to 1 litre with petroleum ether. For preparation of the extraction solution dilute standard solution 1:10 with petroleum ether.

Sample pretreatment: adjust 900 ml water (10 °C) with HCl (12 mol/l) to pH 2 and add 80 g MgSO₄. Add 50 ml of the extraction solution, close the bottle and stir the suspension intensely for 30 min.

Add enough dist. water to separate the organic from the aqueous phase.

Column conditioning: 5 ml petroleum ether


Sample application:

slowly aspirate or force the sample through the column
Elution: wash with 10 ml petroleum ether. Evaporate the combined solutions from sample application and elution to 1 ml at about 75 °C. If necessary, fill up to 1 ml again. (If the hydrocarbon content is high, evaporation to 1 ml may not be necessary.)

Recovery rate: must be > 80 % for *n*-tetracontane.

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Column conditioning: 1 ml *n*-hexane

Sample application:

apply 250 µl waste oil to the column and aspirate or force it into the adsorbent with 2 x 1 ml *n*-hexane

MN Appl. No. 301390

Elution: aspirate or force another 2 x 500 µl *n*-hexane through the column; collect all *n*-hexane fractions and if necessary adjust to a concentration suitable for subsequent analysis by either evaporating *n*-hexane in a stream of nitrogen or by dilution with *n*-hexane

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
Protection from cross contamination

For special applications, which require maximum protection from cross contamination we supply chrome-plated brass valves and stainless steel or PTFE connectors, the application of which is shown below. These special connectors are fitted through the lid; thus the sample only has contact with the inert connector and not with the lid, directly flowing into the receptacle.

Drying attachment

If the eluate has to be evaporated, this can be performed with the so-called drying attachment (**11**, see below). This special lid has a gas connector on one side (**12**), from which the gas is fed simultaneously to the 12 or 24 stations (**13**). Thus 12 or 24 eluates can be evaporated simultaneously by just changing the lid and applying a stream of inert gas, e.g. nitrogen.