

Reagents

You will find here a variety of high quality Avantor DNA Polymerases optimised for specific applications (Taq, glycerol-free Taq, Hot start and High Fidelity enzymes) and ready to use master mixes. Complete your PCR reaction workflow needs with specific Avantor PCR reagents as dNTPs, PCR buffers, PCR clean-up, DNA loading buffer and DNA ladders for electrophoresis gels



Application Guide

	Taq DNA Polymerase	Glycerol-free Taq DNA Polymerase	Red Taq DNA Polymerase	Taq Plus 2X Master Mix	Taq DNA Polymerase Master Mix	Red Taq DNA Polymerase Master Mix	TEMPase Hot Start DNA Polymerase	Glycerol-free TEMPase Hot Start DNA Polymerase	TEMPase Hot Start Master Mix K + C	Blue TEMPase Hot Start Master Mix K + C	GC-rich DNA Target kit	GC TEMPase Master Mix I + II	Multiplex TEMPase Master Mix	Fast HiFi DNA Polymerase	Fast HiFi DNA Polymerase 2X Master Mix
Application	Standard PCR						Hot Start			Special PCR			Hi Fi		
Routine PCR	●		●	●	●	●	●	●	●	●					
High throughput	●	●	●	●	●	●	●	●	●	●					
Automation		●						●							
GC-rich DNA templates							●				●	●		●	●
Multiplex PCR							●						●		
Sequencing														●	●
Genotyping	●	●	●	●	●	●	●	●	●	●			●		
Cloning	●						●							●	●
Mutagenesis														●	●
Freeze-drying		●						●						●	●
Low abundance targets							●	●	●	●	●	●	●		
Forensics													●		
DNA fingerprinting													●		
Colony PCR	●	●	●	●	●	●	●	●	●	●	●	●			
Gene expression							●	●	●	●	●	●			
Microbial detection							●	●	●	●	●	●	●		
NGS applications														●	●
Product example	733-1301	733-1817	733-1323	733-2598	733-2546	733-2542	733-1331	733-2553	733-2581	733-2585	733-2567	733-2561	733-2569	733-2882	733-2886

● Recommended ● Suitable

Technical Guide

	Taq DNA Polymerase	Glycerol-free Taq DNA Polymerase	Red Taq DNA Polymerase	Taq Plus 2X Master Mix	Taq DNA Polymerase Master Mix	Red Taq DNA Polymerase Master Mix	TEMPase Hot Start DNA Polymerase	Glycerol-free TEMPase Hot Start DNA Polymerase	TEMPase Hot Start Master Mix K + C	Blue TEMPase Hot Start Master Mix K + C	GC-rich DNA Target kit	GC TEMPase Master Mix I + II	Multiplex TEMPase Master Mix	Fast HiFi DNA Polymerase	Fast HiFi DNA Polymerase 2X Master Mix
Feature	Standard PCR						Hot Start			Special PCR			Hi Fi		
Direct gel loading						●				●					
Pipetting visualisation			●			●				●					
Proof reading activity														●	●
dUTP incorporation	●	●	●	●	●	●	●	●	●	●	●	●	●		
3'- A overhang	●	●	●	●	●	●	●	●	●	●	●	●	●		
Technical data															
Fidelity versus Taq				1X				1X			< 1X	< 1X	1X		< 60X
Amplicon size				≤ 5 kb				≤ 5 kb				≤ 5 kb			≤ 11.5 kb
Elongation speed				35-100 nt/s				35-100 nt/s				35-100 nt/s			75-100 nt/s
Processivity				60 nt				60 nt				60 nt			
5'-3' exonuclease activity	●	●	●	●	●	●	●	●	●	●	●	●	●		
Performance															
Fidelity	+	+	+	+	+	+	+	+	+	+	+	+	+	+++	+++
Specificity	+	+	+	++	+	+	++	++	++	++	++	++	++	++	++
Sensitivity	+	+	+	+	+	+	++	++	++	++	++	++	++	+	+
Yield	++	++	++	++	++	++	++	++	++	++	++	++	++	++	++

+ High ++ Very High +++ Excellent

TAQ DNA POLYMERASE



VWR® Taq DNA Polymerase is an ultra-pure, thermostable, recombinant DNA polymerase, which provides robust PCR performance in a wide range of PCR applications, without time-consuming optimisation. The enzyme is isolated from *Thermus aquaticus*, and has a molecular weight of approximately 94 kDa. VWR® Taq DNA Polymerase has both a 5' to 3' DNA polymerase and a double strand 5' to 3' exonuclease activity. It leaves an A overhang, which makes the enzyme ideal for TA cloning.

- Ideal choice for routine applications
- High performance, thermostable DNA polymerase
- Optimal for TA cloning

Taq DNA polymerase concentration: 5 Units/μl

10X Key Buffer: Tris-HCl pH 8,5; (NH₄)₂SO₂, 15 mM MgCl₂, 1% Tween® 20

10X Extra Buffer: Tris-HCl pH 8,3; KCl, 15 mM MgCl₂, 1% Triton™ X-100

10X Mg-Free Key Buffer: Tris-HCl pH 8,5; (NH₄)₂SO₂, 1% Tween® 20

10X Mg-Free Extra Buffer: Tris-HCl pH 8,3; KCl, 1% Triton™ X-100

EU = Units

VWR® Taq DNA Polymerase is usually supplied with either or both Key Buffer and Extra Buffer. Key Buffer (NH₄⁺) gives a superior amplification signal (high yield) minimising the need for optimisation of the Mg²⁺ concentration, or the annealing temperature in most primer-template systems. Extra Buffer is a traditional potassium (K⁺) buffer. Extra Buffer promotes high specificity, but careful optimisation of primer annealing temperatures and Mg²⁺ concentrations may be required.

Description	Taq DNA Polymerase 5U/μl 250 units	Taq DNA Polymerase 5U/μl 500 units	Taq DNA Polymerase 5U/μl 1000 units	Taq DNA Polymerase 5U/μl 2.500 units	Taq DNA Polymerase 5U/μl 5.000 units	Taq DNA Polymerase 5U/μl 10000 units	Taq DNA Polymerase 5U/μl 1000 units	Taq DNA Polymerase 5U/μl 2.500 units	Taq DNA Polymerase 5U/μl 10.000 units	Taq DNA Polymerase 5U/μl 1000 units	Taq DNA Polymerase 5U/μl 10.000 units	
Amplicon size	≤ 5kb											
Application	Standard amplification											
Buffer composition	10x Key Buffer, 10x Extra Buffer and 25 mM MgCl ₂						10x Key Buffer Mg-free and 25 mM MgCl ₂			10x Extra Buffer and 25 mM MgCl ₂ Detergent free	10x Key Buffer and 25 mM MgCl ₂ Detergent free	10x Key Buffer, 10x Extra Buffer and 25 mM MgCl ₂ Detergent free
Concentration	5 U/μl											
Exonuclease activity	5' → 3'											
Features							Mg-free			Detergent free		
Format	Enzyme with/without buffer											
High fidelity	No											
Hot start												
Processivity	Standard											
Product overhang	3' A											
Size	1 x 0.05 ml + 3 x 1.5 ml	1 x 0.1 ml + 3 x 1.5 ml	2 x 0.1 ml + 6 x 1.5 ml	5 x 0.1 ml + 15 x 1.5 ml	10 x 0.1 ml + 9 x 5 ml	3 x 0.667 ml + 18 x 5 ml	2 x 0.1 ml + 4 x 1.5 ml	5 x 0.1 ml + 10 x 1.5 ml	3 x 0.667 ml + 12 x 5 ml	2 x 0.1 ml + 4 x 1.5 ml	3 x 0.667 ml + 18 x 5 ml	
Storage conditions	20 °C (long term), + 4 °C (6 months)											
Units	250	500	1000	2500	5000	10000	1000	2500	10000	1000	10000	

Description	Pk	Cat. No.
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	250 EU	733-1300
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	500 EU	733-1301
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	1.000 EU	733-1302
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	2.500 EU	733-1819
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	5.000 EU	733-1820
Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	10.000 EU	733-1303
Taq DNA Polymerase, 10X MgCl ₂ -free Key Buffer, 25 mM MgCl ₂	1.000 EU	733-1312
Taq DNA Polymerase, 10X MgCl ₂ -free Key Buffer, 25 mM MgCl ₂	2.500 EU	733-1313
Taq DNA Polymerase, 10X MgCl ₂ -free Key Buffer, 25 mM MgCl ₂	10.000 EU	733-2009
Taq DNA Polymerase, 10X MgCl ₂ -free Extra Buffer, 25 mM MgCl ₂	1.000 EU	733-1305
Taq DNA Polymerase, 10X Tween-free Key Buffer (15 mM MgCl ₂), 10X Triton-free Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	1.000 EU	733-1307
Taq DNA Polymerase, 10X Tween-free Key Buffer (15 mM MgCl ₂), 10X Triton-free Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	10.000 EU	733-1823

TAQ DNA POLYMERASE MASTER MIX



VWR® Taq DNA Polymerase Master Mix is a ready to use 1,1X or 2X reaction mix. Simply add primers, template and water to carry out primer extensions and other molecular biology applications.

Tests = Reactions

Description	Taq DNA Polymerase 1,1X Master M 1,5 mM MgCl ₂ 500 reactions	Taq DNA Polymerase 1,1X Master Mix 1,5 mM MgCl ₂ 2500 reactions	Taq DNA Polymerase 1,1X Master Mix 2 mM MgCl ₂ 2500 reactions	Taq DNA Polymerase 2X Master Mix 1,5 mM MgCl ₂ 500 reactions	Taq DNA Polymerase 2X Master Mix 1,5 mM MgCl ₂ 2500 reactions	Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 500 reactions	Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 2500 reactions
Amplicon size	≤ 5kb						
Application	Standard amplification						
Concentration	1.1x			2x			
Exonuclease activity	5' →3'						
Format	Mastermix incl. dNTPs						
High fidelity	No						
Hot start	No						
Processivity	Standard						
Product overhang	3' A						
Size	15 x 1.5 ml	75 x 1.5 ml		10 x 1.25 ml	50 x 1.25 ml	10 x 1.25 ml	50 x 1.25 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)						
Units	500 r Mix 500	2500 r Mix 2500		500 r Mix 500	2500 r Mix 2500	500 r Mix 500	2500 r Mix 2500

Description	Pk	Cat. No.
Taq DNA Polymerase 1,1X Master Mix, 1,5 mM MgCl ₂	500 Tests	733-2540
Taq DNA Polymerase 1,1X Master Mix, 1,5 mM MgCl ₂	2.500 Tests	733-1314
Taq DNA Polymerase 1,1X Master Mix, 2,0 mM MgCl ₂	2.500 Tests	733-1315
Taq DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	500 Tests	733-2542
Taq DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	2.500 Tests	733-1316
Taq DNA Polymerase 2X Master Mix, 2,0 mM MgCl ₂	500 Tests	733-2543
Taq DNA Polymerase 2X Master Mix, 2,0 mM MgCl ₂	2.500 Tests	733-1317

VWR® TAQ DNA POLYMERASE GLYCEROL FREE



VWR Taq DNA Polymerase Glycerol Free 5U/μl is a thermostable recombinant DNA polymerase, which exhibits very high activity in primer extension and other molecular biology applications. The enzyme is isolated from *Thermus aquaticus* and has a molecular weight of approximately 94 kDa.

VWR Taq DNA Polymerase Glycerol Free has both a 5'→3' DNA polymerase and a 5'→3' exonuclease activity. The enzyme lacks a 3'→5' exonuclease activity (no proofreading ability). VWR Taq DNA Polymerase Glycerol Free leaves an A' overhang, which makes the enzyme ideal for TA cloning.

Description	Taq DNA Polymerase Glycerol Free 5U/μl 100.000 units	Taq DNA Polymerase Glycerol Free 5U/μl 375.000 units	Taq DNA Polymerase Glycerol Free 5U/μl 2.000.000 units	Taq DNA Polymerase Glycerol Free 5U/μl 5.000 units	Taq DNA Polymerase Glycerol Free 5U/μl 1.000 units
Amplicon size	≤ 5kb				
Application	Standard amplification, Freeze-drying				
Buffer composition	No buffer			25 mM MgCl ₂	10x Key Buffer, 10x Extra Buffer and 25 mM MgCl ₂
Concentration	5 U/μl				
Exonuclease activity	5' →3'				
Features	Glycerol free				
Format	Enzyme with/without buffer				
High fidelity	No				
Hot start	No				
Processivity	Standard				
Product overhang	3' A				
Size	20 x 1 ml	1 x 75 ml	Bulk	10 x 0.1 ml	2 x 0.1 ml + 6 x 1.5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)				
Units	100000	375000	2000000	5000	1000

Description	Pk	Cat. No.
Taq DNA Polymerase Glycerol Free 5U/μl 100.000 units	20ml	733-2406
Taq DNA Polymerase Glycerol Free 5U/μl 375.000 units	75 ml	733-3843
Taq DNA Polymerase, glycerol-free, 200 000 units	200.000 EU	733-2038
Taq DNA Polymerase, glycerol-free, 25 mM MgCl ₂	5.000 EU	733-1999
Taq DNA Polymerase, glycerol-free, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	1.000 EU	733-1817

VWR® RED TAQ DNA POLYMERASE

VWR® Red Taq DNA Polymerase is a blend of Taq DNA polymerase combined with an inert red dye. The dye enables quick visual recognition of reactions to which enzyme has been added, as well as confirmation of complete mixing.



- Ideal choice for routine applications
- High performance, thermostable DNA polymerase
- Optimal for TA cloning

Taq DNA polymerase concentration: 5 Units/μl

10X Key Buffer: Tris-HCl pH 8,5; (NH₄)₂SO₂, 15 mM MgCl₂, 1% Tween® 20

EU = Units

VWR® Taq DNA Polymerase is usually supplied with Key Buffer. Key Buffer (NH₄⁺) gives a superior amplification signal (high yield) minimising the need for optimisation of the Mg²⁺ concentration, or the annealing temperature in most primer-template systems.

Description	"Taq RED DNA Polymerase 5U/μl 500 units"	"Taq RED DNA Polymerase 5U/μl 1.000 units"	"Taq RED DNA Polymerase 5U/μl 2.500 units"	"Taq RED DNA Polymerase 5U/μl 10.000 units"
Amplicon size	≤ 5kb			
Application	Standard amplification			
Buffer composition	10x Key Buffer, 10x Extra Buffer and 25 mM MgCl ₂			"10x Key Buffer, 10x Extra Buffer and M25 mM gCl ₂ "
Concentration	5 U/μl			
Exonuclease activity	5' →3'			
Features	Visible dye			
Format	Enzyme with/without buffer			
High fidelity	No			
Hot start				
Processivity	Standard			
Product overhang	3' A			
Size	1 x 0.1 ml + 3 x 1.5 ml	2 x 0.1 ml + 6 x 1.5 ml	5 x 0.1 ml + 15 x 1.5 ml	3 x 0.667 ml + 18 x 5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)			
Units	500	1000	2500	10000

Description	Pk	Cat. No.
Red Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	500 EU	733-2408
Red Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	1.000 EU	733-2409
Red Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	2.500 EU	733-1323
Red Taq DNA Polymerase, 10X Key Buffer (15 mM MgCl ₂), 10X Extra Buffer (15 mM MgCl ₂), 25 mM MgCl ₂	10.000 EU	733-1834



RED *Taq* DNA POLYMERASE MASTER MIX

VWR® Red *Taq* DNA Polymerase Master Mix, which also contains an inert red dye, can be directly loaded onto an agarose gel without addition of electrophoresis loading buffers.

Tests = Reactions

Description	Red Taq DNA Polymerase 1.1X Master Mix 1.5 mM MgCl ₂ 500 reactions	Red Taq DNA Polymerase 1.1X Master Mix 1.5 mM MgCl ₂ 2500 reactions	Red Taq DNA Polymerase 1.1X Master Mix 2 mM MgCl ₂ 500 reactions	Red Taq DNA Polymerase 1.1X Master Mix 2 mM MgCl ₂ 2500 reactions	Red Taq DNA Polymerase 2X Master Mix 1.5 mM MgCl ₂ 500 reactions	Red Taq DNA Polymerase 2X Master Mix 1.5 mM MgCl ₂ 2500 reactions	Red Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 5000 reactions	Red Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 10.000 reactions	Red Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 20 000 reactions	Red Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 500 reactions	Red Taq DNA Polymerase 2X Master Mix 2 mM MgCl ₂ 2500 reactions
Amplicon size	≤ 5kb										
Application	Standard amplification										
Concentration	1.1x					2x					
Exonuclease activity	5' → 3'										
Features	Direct loading										
Format	Mastermix incl. dNTPs										
High fidelity	No										
Hot start	Standard										
Processivity	3' A										
Product overhang	3' A										
Size	15 x 1.5 ml	75 x 1.5 ml	15 x 1.5 ml	75 x 1.5 ml	10 x 1.25 ml	50 x 1.25 ml	25 x 5 ml	28 x 9 ml	1 x 500 ml	10 x 1.25 ml	50 x 1.25 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)										
Units	500 r Mix 500	2500 r Mix 2500	500 r Mix 500	2500 r Mix 2500	500 r Mix 500	2500 r Mix 2500	5000 r Mix 5000	10000 r Mix 10000	20000 r Mix 20 000	500 r Mix 500	2500 r Mix 2500

Description	Pk	Cat. No.
Red <i>Taq</i> DNA Polymerase 1,1X Master Mix, 1,5 mM MgCl ₂	500 Tests	733-2544
Red <i>Taq</i> DNA Polymerase 1,1X Master Mix, 1,5 mM MgCl ₂	2.500 Tests	733-1318
Red <i>Taq</i> DNA Polymerase 1,1X Master Mix, 2,0 mM MgCl ₂	500 Tests	733-2545
Red <i>Taq</i> DNA Polymerase 1,1X Master Mix, 2,0 mM MgCl ₂	2.500 Tests	733-1319
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	500 Tests	733-2546
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	2.500 Tests	733-1320
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	5.000 Tests	733-2130
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	10.000 Tests	733-2131
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 1,5 mM MgCl ₂	20.000 Tests	733-2132
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 2,0 mM MgCl ₂	500 Tests	733-2547
Red <i>Taq</i> DNA Polymerase 2X Master Mix, 2,0 mM MgCl ₂	2.500 Tests	733-1321



Taq Plus 2x Master Mix

TAQ PLUS 2X MASTER MIXES

VWR® *Taq* Plus is an optimised format of *Taq* DNA polymerase master mix and, therefore, is a good alternative to *Taq* DNA polymerase and *Taq* DNA polymerase master mix.

- Ready to use 2X master mix for greater convenience
- Significantly reduced set-up time
- Increased specificity compared to *Taq* DNA polymerase master mix
- Reduced risk of contamination

Ensures increased specificity and improved PCR performance, even on difficult templates. These reagents are suitable for standard PCR applications, screening and high throughput testing.

This PCR master mix is composed of *Taq* DNA polymerase and an optimised buffer system, including dNTPs and magnesium chloride (1,5 mM).

Description	Taq Plus 2x Master Mix 1.5 mMNgCl2 100 reactions	Taq Plus 2x Master Mix 1.5 mMNgCl2 500 reactions	Taq Plus 2x Master Mix 1.5 mMNgCl2 2500 reactions	Taq Plus 2x Master Mix 1.5 mMNgCl2 5000 reactions
Amplicon size	≤ 4kb			
Application	Standard amplification			
Concentration	2x			
Exonuclease activity	5' →3'			
Format	Mastermix incl. dNTPs			
High fidelity	No			
Hot start	Standard			
Processivity	3' A			
Product overhang	3' A			
Size	2 x 1.25 ml	10 x 1.25 ml	50 x 1.25 ml	25 x 5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)			
Units	100	500	2500	5000

Description	Pk	Cat. No.
Taq Plus 2X master mix, 2x1,25 ml	100 Tests	733-2597
Taq Plus 2X master mix, 10x1,25 ml	500 Tests	733-2598
Taq Plus 2X master mix, 50x1,25 ml	2.500 Tests	733-2599
Taq Plus 2X master mix, 25x5 ml	5.000 Tests	733-2600



TEMPASE HOT START DNA POLYMERASE

VWR® TEMPase Hot Start DNA Polymerases are highly stable polymerases, featuring higher specificity, superior sensitivity and greater yields compared to standard DNA polymerases. These features make them well suited for the detection of low abundance targets. Other uses include screening, amplification of GC-rich sequences, multiplex PCR, direct PCR and qPCR. A glycerol-free TEMPase Hot Start DNA Polymerase is also available for automation and freeze drying.

The GC-Rich Template kit is specifically designed for difficult GC-rich sequences. Combined with TEMPase, GC buffers I and II promote excellent amplification. The kit is designed for initial testing before using one of the GC-TEMPase 2X Master Mixes.

VWR® TEMPase DNA polymerases generally include two different buffers, Key Buffer and Combination Buffer, which are each suited to different PCR requirements. Key Buffer (NH⁴⁺) gives a superior amplification signal (high yield), and minimises the need for optimisation of the Mg²⁺ concentration or the annealing temperature in most primer-template systems. Combination Buffer is a mixture of K⁺ and NH⁴⁺. It combines high specificity with good product yield and high tolerance to optimisation of primer annealing temperatures and Mg²⁺ concentrations, due to its balanced ammonium-potassium formulation. Each buffer contains 15 mM MgCl₂ (1,5 mM in final volume). Additional MgCl₂ for easy optimisation is included in a separate vial.

Description	TEMPase HotStart DNA Polymerase 5U/μl 500 units	TEMPase HotStart DNA Polymerase 5U/μl 2500 units	TEMPase HotStart DNA Polymerase 5U/μl 10.000 units	TEMPase HotStart DNA Polymerase Glycerol Free 5U/μl 500 units	TEMPase HotStart DNA Polymerase Glycerol Free 5U/μl 2500 units	
Amplicon size	≤ 5kb					
Application	Standard amplification, Low-quality template, Multiplex			Standard amplification, Freeze-drying		
Buffer composition	10x Key Buffer, 10x Combination Buffer and 25 mM MgCl ₂			25 mM MgCl ₂	10x Key Buffer and 25 mM MgCl ₂	10x Combination Buffer and 25 mM MgCl ₂
Concentration	5 U/μl					
Exonuclease activity	5' →3'					
Features	Glycerol free					
Format	Enzyme with/without buffer					
High fidelity	No					
Hot start	Chemically modified					
Processivity	Standard					
Product overhang	3' A					
Size	1 x 0.1 ml + 3 x 1.5 ml	5 x 0.1 ml + 15 x 1.5 ml	3 x 0.667 ml + 18 x 5 ml	1 x 0.1 ml + 1 x 1.5 ml	5 x 0.1 ml + 5 x 1.5 ml	5 x 0.1 ml + 10 x 1.5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)					
Units	500	2500	10000	500	2500	

Description	Pk	Cat. No.
TEMPase Hot Start DNA Polymerase, 5 U/μl, with 10X Key Buffer, 10X Combination Buffer and MgCl ₂	500 EU	733-1331
TEMPase Hot Start DNA Polymerase, 5 U/μl, with 10X Key Buffer, 10X Combination Buffer and MgCl ₂	2.500 EU	733-1333
TEMPase Hot Start DNA Polymerase, 5 U/μl, with 10X Key Buffer, 10X Combination Buffer and MgCl ₂	10.000 EU	733-1838
TEMPase Hot Start DNA Polymerase, glycerol-free, 5 U/μl, without buffers	500 Tests	733-2552
TEMPase Hot Start DNA Polymerase, glycerol-free, 5 U/μl, without buffers	2.500 Tests	733-2553
TEMPase Hot Start DNA Polymerase, glycerol-free, 5 U/μl, with 10X Key Buffer and MgCl ₂	2.500 Tests	733-2556
TEMPase Hot Start DNA Polymerase, glycerol-free, 5 U/μl, with 10X Combination Buffer and MgCl ₂	2.500 Tests	733-2559



TEMPASE HOT START 2x MASTER MIX

TEMPase Hot Start DNA Polymerase Master Mix and Blue TEMPase Master Mix offer easy reaction assembly at room temperature, reduced set-up time and fewer handling steps, which lead to increased reproducibility. As a consequence TEMPase Hot Start DNA Polymerase Master Mix is highly suited to standard tests.

The blue loading dye in Blue TEMPase Hot Start DNA Polymerase Master Mix facilitates direct gel loading, and eliminates the need for separate loading dye - no need for time-consuming sample preparation before electrophoresis.

Multiplex 2x Master Mix is composed of TEMPase Hot Start DNA Polymerase and a specialised buffer system designed for multiplex PCR.

TEMPase Hot Start DNA Polymerase Master Mix and Blue TEMPase Master Mix are available in two variations, either based on Key Buffer (Master Mix K) or Combination Buffer (Master Mix C) to suit different PCR requirements. Additional MgCl₂ is included in the kit to enable optimisation.

Description	TEMPase Hot Start DNA Polymerase 2X Master Mix A 1.5 mM MgCl ₂ 2500 reactions	TEMPase Hot Start DNA Polymerase 2X Master Mix C 1.5 mM MgCl ₂ 2500 reactions	Blue TEMPase Hot Start DNA Polymerase 2X Master Mix A 1.5 mM MgCl ₂ 2500 reactions	Blue TEMPase Hot Start DNA Polymerase 2X Master Mix C 1.5 mM MgCl ₂ 2500 reactions	Multiplex TEMPase 2X Master Mix 1.5 mM MgCl ₂ 2500 reactions
Amplicon size	≤ 5 kb				
Application	Standard amplification, Low-quality template, Multiplex				Standard amplification, Multiplex
Buffer composition	Based on Ammonium buffer	Based on Combination buffer	Based on Ammonium buffer	Based on Combination buffer	25 mM MgCl ₂
Concentration	2x				
Exonuclease activity	5' → 3'				
Features	Direct loading				
Format	Mastermix incl. dNTPs				
High fidelity	No				
Hot start	Chemically modified				
Processivity	Standard				
Product overhang	3' A				
Size	50 x 1.25 ml				50 x 1.25 ml + 3 x 1.5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)				
Units	2500 r Mix 2500				

Description	Pk	Cat. No.
TEMPase Hot Start 2X Master Mix, with Master Mix K, 1,5 mM MgCl ₂	2.500 Tests	733-2582
TEMPase Hot Start 2X Master Mix, with Master Mix C, 1,5 mM MgCl ₂	2.500 Tests	733-1840
Blue TEMPase Hot Start 2X Master Mix, with Blue Master Mix K, 1,5 mM MgCl ₂	2.500 Tests	733-2585
Blue TEMPase Hot Start 2X Master Mix, with Blue Master Mix C, 1,5 mM MgCl ₂	2.500 Tests	733-1841
Multiplex TEMPase Hot Start 2X Master Mix, 1,5 mM MgCl ₂ with separate vial of MgCl ₂	2.500 Tests	733-2569



Fast HiFi DNA polymerase 2 U/μl

FAST HIFI DNA POLYMERASE, 2 U/ML

Fast HiFi DNA Polymerase 2 U/μl is a proofreading DNA polymerase displaying the following features; high fidelity >60X Taq DNA Polymerase, ability to amplify problematic DNA targets, such as those with low to high GC content and ability to perform amplification on long DNA targets. It is recommended for applications, which require extremely high fidelity such as cloning/sub-cloning, NGS applications and mutagenesis.

- High fidelity: >60X Taq Fidelity
- 3'→5' proofreading exonuclease activity
- Processes up to 18 kb from complex templates (eg. Human gDNA)
- Good coverage on DNA templates with low to high GC content
- Fast extension rate: 10 s/kb

Fast HiFi DNA Polymerase 2 U/μl exhibits both 5'→3' DNA polymerase activity and 3'→5' proof reading exonuclease activity enabling this polymerase to correct base pair mismatches. It is supported with 5X Fast HiFi Buffer allowing robust amplification on DNA targets with low to high GC content and long DNA targets.

For difficult amplicons, such as GC-rich DNA samples, those with complex secondary structures or long amplicons, the addition of 1 to 2 M Betaine enhancer solution is recommended.

For more convenient handling Fast HiFi DNA Polymerase 2 U/μl is also available as a 2X master mix.

Sample types: gDNA, plasmid DNA, bacterial DNA, λDNA and other DNA templates

Description	Fast HiFi DNA Polymerase 2 U/μl 100 units	HiFi DNA Polymerase 2 U/μl 500 units	Fast HiFi DNA Polymerase 2 U/μl 1000 units	Fast HiFi DNA Polymerase 2 U/μl 2500 units
Amplicon size	18 kb for gDNA			
Application	Standard amplification, Cloning (high fidelity), Long PCR, NGS library amplification			
Buffer composition	10X Fast HiFi Buffer + 25 mM MgCl			
Concentration	2 U/μl			
Exonuclease activity	5' →3' 3' →5'			
Features	Proof ready activity, Long template			
Format	Enzyme with/without buffer			
High fidelity	Yes			
Hot start	No			
Processivity	Fast			
Product overhang	No			
Size	1x0,05 ml + 2x1,5 ml	1x0,25 ml + 3x1,5 ml	2x0,25 ml + 6x1,5 ml	5x0,25 ml + 14x1,5 ml
Storage conditions	20 °C (long-term), +4 °C (6 months)			
Units	100	500	1000	2500

Description	Pk	Cat. No.
Fast HiFi DNA Polymerase 2 U/μl, 100 U	100 EU	733-2880
Fast HiFi DNA Polymerase 2 U/μl, 500 U	500 EU	733-2881
Fast HiFi DNA Polymerase 2 U/μl, 1000 U	1.000 EU	733-2882
Fast HiFi DNA Polymerase 2 U/μl, 2500 U	2.500 EU	733-2883



HiFi DNA polymerase 2x master mix

FAST HiFi DNA POLYMERASE, 2X MASTER MIX, VWR®

Fast HiFi DNA Polymerase 2 U/μl is a proofreading DNA polymerase displaying the following features; high fidelity >60X Taq DNA Polymerase, ability to amplify problematic DNA targets, such as those with low to high GC content and ability to perform amplification on long DNA targets. It is recommended for applications, which require extremely high fidelity such as cloning/sub-cloning, NGS applications and mutagenesis.

- All-in-one 2x master mix for great convenience
- High fidelity: >60X Taq Fidelity
- 3'→5' proofreading exonuclease activity
- Processes up to 11 kb from complex templates (eg. Human gDNA)
- Good coverage on DNA templates with low to high GC content
- Fast extension rate: 10 s/kb

Fast HiFi DNA Polymerase 2x master mix is a ready to use 2x reaction mix composed of Fast HiFi DNA Polymerase and an optimised buffer system including dNTPs and Magnesium chloride, allowing robust amplification on DNA target with low to high GC content and long DNA targets. Fast HiFi DNA Polymerase exhibits both 5'→3' DNA polymerase activity and 3'→5' proofreading exonuclease activity enabling this polymerase to correct base pair mismatches.

For difficult amplicons, such as GC-rich DNA samples, those with complex secondary structures or long amplicons, the addition of 1 to 2 M Betaine enhancer solution is recommended.

Fast HiFi DNA Polymerase 2x master mix is recommended for applications, which require extremely high fidelity such as cloning/sub-cloning, NGS applications and mutagenesis.

Sample types: gDNA, plasmid DNA, bacterial DNA, λDNA and other DNA templates

Description	Fast HiFi DNA Polymerase 2x Master Mix 100 reactions	Fast HiFi DNA Polymerase 2x Master Mix 500 reactions	Fast HiFi DNA Polymerase 2x Master Mix 2500 reactions	Fast HiFi DNA Polymerase 2x Master Mix 5000 reactions
Amplicon size	11 kb for gDNA			
Application	Standard amplification, Cloning (high fidelity), Long PCR, NGS library amplification			
Concentration	2x			
Exonuclease activity	5' →3' 3' →5'			
Features	Proof ready activity, Long template			
Format	Mastermix incl. dNTPs			
High fidelity	Yes			
Hot start	No			
Processivity	Fast			
Product overhang	No			
Size	2 x 1.25 ml	10 x 1.25 ml	50 x 1.25 ml	25 x 5 ml
Storage conditions	20 °C (long term), + 4 °C (6 months)			
Units	100 r Mix 100	500 r Mix 500	2500 r Mix 2500	5000 r Mix 5000

Description	Pk	Cat. No.
Fast HiFi DNA Polymerase, 2X master mix, 100 reactions	100 Tests	733-2884
Fast HiFi DNA Polymerase, 2X master mix, 500 reactions	500 Tests	733-2885
Fast HiFi DNA Polymerase, 2X master mix, 1000 reactions	2.500 Tests	733-2886
Fast HiFi DNA Polymerase, 2X master mix, 2500 reactions	5.000 Tests	733-2887



ExoCleanUp

ONE-STEP PCR CLEAN-UP, EXOCLEANUP FAST

VWR ExoCleanUp FAST PCR reagent is a one-step PCR clean-up reagent for optimal sequencing results, consisting of a balanced combination of a heat-labile exonuclease I (HL-Exol) and a recombinant shrimp alkaline phosphatase (rSAP).

- Designed to clean-up PCR products in 5 minutes
- No need for spin columns or magnetic beads
- Treatment improves downstream applications, such as DNA sequencing and SNP analysis

Treatment of amplified PCR products with this reagent helps to remove residual primers and single-stranded DNA, and inactivates excess dNTPs by dephosphorylation. After enzymatic treatment at 37 °C for a minimum of 2 minutes, this reagent is completely inactivated by heating at 80 °C for a minimum of 3 minutes.

Description	ExoCleanUp FAST PCR clean-up reagent			
Method/Format	Enzymatic clean-up			
Prep Size	Scalable			
Refine target molecule	Amplified/modified DNA			
Sample size	1 × 0,2 ml	1 × 1 ml	10 × 1 ml	4 × 1 ml
Sample type	PCR and other enzymatic reactions			
Target molecule	DNA			

Description	Pk	Cat. No.
ExoCleanUp FAST PCR clean-up reagent	100 Tests	733-2592
ExoCleanUp FAST PCR clean-up reagent	500 Tests	733-2593
ExoCleanUp FAST PCR clean-up reagent	5.000 Tests	733-2594
ExoCleanUp FAST PCR clean-up reagent	2.000 Tests	733-2849



Fast extract DNA solution

FAST EXTRACT DNA SOLUTION, VWR®

The fast extract DNA solution provides rapid and efficient extraction of PCR-ready DNA from mammalian tissues. DNA is ready in eight minutes. The one-reagent protocol is divided into two simple heating steps, which can be directly followed by PCR analysis, such as screening and genotyping.

- One-reagent set-up
- Rapid eight-minute protocol
- PCR-ready DNA
- DNA extracts from mammalian tissues and bacteria
- Non-toxic reagents
- Automation-friendly

DNA can be extracted from many different sample types e.g. mouse tails or ears, saliva, bacteria and mammalian tissues. The one-reagent DNA extraction set-up is easily scaled and can be conducted by robotic automation platforms. Depending on the sample size, the DNA extraction can be performed in PCR tubes or in 1,5 ml tubes, using either a thermocycler or a heating block.

Prep size: 100 µl

Method/format: DNA extraction

Description	Pk	Cat. No.
VWR® fast extract DNA solution, 8 min, 100 reactions	100 Tests	733-2876
VWR® fast extract DNA solution, 8 min, 500 reactions	500 Tests	733-2877



Fast extract genotyping PCR kit

FAST EXTRACT GENOTYPING PCR KIT, VWR®

The fast extract genotyping PCR kit consist of fast extraction DNA solution and VWR® Red Taq DNA polymerase 2X Master mix. This kit is ideal for genotyping of DNA extracted fom mammalian tissues, saliva or bacteria and provides PCR-ready DNA in eight minutes and genotyping results in less than one and half an hour.

- Kit for genotyping
- Reagents for DNA extraction and PCR
- Rapid eight-minute extraction protocol
- PCR-ready DNA
- DNA extracts from mammalian tissues
- Red loading dye within PCR master mix for direct gel loading

The Fast Extract DNA Solution provides rapid and efficient extraction of PCR-ready DNA from mammalian tissues. DNA is ready in eight minutes. The one-reagent protocol is divided into two simple heating steps, which can be directly followed by PCR analysis using Red Taq DNA Polymerase 2X Master Mix. DNA extraction set-up is easily scaled and can be conducted by robotic automation platforms. Depending on the sample size, the DNA extraction can be performed in PCR tubes or in 1,5 ml tubes, using either a thermocycler or a heating block.

Red Taq DNA polymerase 2X Master mix is a ready to use 2X reaction mix. Simply add primers, the extracted DNA to the PCR reaction mix to successfully carry out PCR. The red dye and stabiliser within the master mix formulation enables the user to load amplified PCR samples directly to the DNA gel. The red dye also provides visualisation of pipetting and mixing.

Prep size: 100 + 25 µl

Method/format: Genotyping - DNA extraction + PCR

Description	Pk	Cat. No.
VWR® fast extract genotyping PCR kit, 100 reactions	100 Tests	733-2878
VWR® fast extract genotyping PCR kit, 500 reactions	500 Tests	733-2879

WATER, PCR GRADE

Ultrapure water, free of endonuclease-, nicking-, and exonuclease activity, free of human DNA.

- Convenient 5 ml aliquots
- Tested for contaminating activities

Storage at room temperature. Long term storage at -20 °C. Product expiry at -20 °C is stated on the label.

Description	Pack type	Pk Info	Pk	Cat. No.
Water, PCR grade	Plastic tube	6x5 ml	1 KIT	733-2573



dNTP mix, 10 mM

DNTPS

Ready to use molecular biology grade dNTP mixes and dNTP sets. The dNTP mix is designed to save hands-on time for researchers, and reduce the possibility of contamination by reducing pipetting. dNTP solutions are also available in sets of four individual dNTPs, each 100 mM. Both are convenient for use in DNA polymerisation reactions, DNA labelling and sequencing processes.

- Available as pre-mixed 10 or 25 mM solutions, or as sets of individual 100 mM dNTP solutions
- Both pre-mixed and sets have been functionally tested in PCR
- Purity >99% by HPLC
- Supplied in solution at pH 7,3 to 7,5

Description	Pk	Cat. No.
dNTP mix, 10 mM of each dA, dC, dG, and dT, 2x500 µl	1.000 µl	733-1363
dNTP set, separate vials of dA, dC, dG, dT, each 100 mM, 4x250 µl	1 SET	733-1364
dNTP mix, 25 mM of each dA, dC, dG, and dT, 2x1 ml	2.000 µl	733-1854
dNTP set, separate vials of dA, dC, dG, dT, each 100 mM, 16x250 µl	1 SET	733-1855



10X Key Buffer

PCR BUFFERS

An optimal buffer system is critical for the performance of successful PCR. VWR® Taq DNA Polymerase kits generally include two different buffers, Key Buffer and Extra Buffer, which are suited for different PCR needs. All buffers contain Tris and 15 mM MgCl₂ (1,5 mM MgCl₂ final concentration). Additional MgCl₂ for easy optimisation is included in a separate vial.

Key Buffer: Key Buffer (NH⁴⁺) gives a superior amplification signal (high yield) and minimises the need for optimisation of the Mg²⁺ concentration or the annealing temperature in most primer-template systems.

Extra Buffer: Extra Buffer is the traditional potassium (K⁺) buffer. Extra Buffer promotes high specificity. Careful optimisation of primer annealing temperatures and Mg²⁺ concentrations may be required.

Combination Buffer: Combination Buffer is a proprietary mixture of K⁺ and NH⁴⁺. It combines high specificity with good product yield and high tolerance to optimisation of primer annealing temperatures and Mg²⁺ concentrations due to its balanced ammonium-potassium formulation.

Buffers for GC-rich templates: Combined with VWR® TEMPase the GC buffers promote excellent amplification results with targets of varying degrees of GC content. This confirms that the choice of buffer is crucial for any successful amplification.

Description	Key Buffer	Extra Buffer	Combination Buffer
Application		Standard amplification	
Concentration		10x	
Size		3 x 1.5 ml	
Storage conditions		20 °C (long term), + 4 °C (6 months)	

Description	Pk	Cat. No.
10X Key Buffer, 15 mM MgCl ₂	1 SET	733-1349
10X Extra Buffer, 15 mM MgCl ₂	1 KIT	733-2303
10X Combination Buffer, 15 mM MgCl ₂	1 SET	733-1352



5X Betaine enhancer solution

BETAINE ENHANCER, 5 M (5X)

Betaine enhancer is especially effective when used with high GC-rich regions or templates with a high degree of secondary structures. It has a decreasing effect on the primer melting temperature.

Description	Pk	Cat. No.
5X Betaine enhancer solution, 5 M, 5x1 ml	5 ml	733-1361
5X Betaine Enhancer solution, 5 M, 100x10 ml	1 KIT	733-2450
5X Betaine Enhancer solution, 5 M, 10x100 ml	1 KIT	733-2451



DNA loading dyes

DNA LOADING BUFFERS

These DNA loading dyes are used to load DNA samples to agarose or SDS DNA gels for gel electrophoresis.


Suitable for TAE, TBE, SB and LB electrophoresis buffers.

EC Regulation No. 1907/2006 (REACH).

Each kit contains 5x1 ml.

Store at -20 °C for *in vitro* use only.

Description	Colour	Pk	Cat. No.
Loading dye 5X	Red	1 KIT	733-2574
Loading dye 5X	Blue	1 KIT	733-2575
Loading dye 5X	Orange	1 KIT	733-2576



DNA ladder 100 BP

DNA LADDERS FOR AGAROSE AND POLYACRILAMIDE GELS

VWR DNA ladders are supplied in a loading buffer that is ready to use on agarose and polyacrylamide DNA gels. The ladders are suitable with TBE, TAE, SB and LB electrophoresis systems.

- Supplied in loading buffer, ready to use directly on the gel
- Mass-calibrated bands for DNA quantification
- PCR DNA ladder has a 1000 bp band that is extra bright to serve as reference point


High range DNA ladder molecular range: 200 to 12000 bp; mass-calibrated bands from 6 to 40 ng/band for DNA quantification.

Low range DNA ladder molecular range: 100 to 1000 bp; mass-calibrated bands from 8 to 40 ng/band for DNA quantification.

PCR DNA ladder molecular range: 100 to 3000 bp; mass-calibrated bands from 10 to 30 ng/band for DNA quantification.

Each DNA ladder is supplied in 0,5 ml packs sufficient for 250 lanes.

Description	Pk	Cat. No.
High range DNA ladder, 200 to 12000 bp, sufficient for 250 lanes (loading 2 µl per gel), 1x0,5 ml	1 KIT	733-2577
Low range DNA ladder, 100 to 1000 bp, sufficient for 250 lanes (loading 2 µl per gel), 1x0,5 ml	1 KIT	733-2578
PCR ladder, 100 to 3000 bp, sufficient for 250 lanes (loading 2 µl per gel), 1x0,5 ml	1 KIT	733-2579



Low DNA ladder

DNA LADDERS FOR AGAROSE AND SDS GELS

VWR DNA ladders are convenient dsDNA ladders supplied in a loading buffer, which are ready to use on agarose and SDS DNA gels.

- Ready to use loading buffer
- For direct loading and easy visualisation
- Suitable with TBE and TAE electrophoresis systems
- Clear, distinct bands
- Blue dye front running at 100 to 300 bp at 0,5 to 1,5% agarose

Supplied in 0,5 ml packs for 100 lanes.

Description	Pk	Cat. No.
Mini DNA ladder, 100 – 500 bp	100 Tests	733-2601
Low range DNA ladder, 100 – 1000 bp	100 Tests	733-2602
PCR ladder, 100 – 3000 bp	100 Tests	733-2603

From sample to sequence

DISCOVER
MORE AT
VWR.COM

From manual to automated processes


Increase your throughput and gain more insight while saving costs and time.



PCR consumables

Good quality at reasonable price?
Easy to find at **vwr.com**. Will it fit
your application or PCR instrument?
Use these selector tools!

Quick search
by selector



Click here to select **PCR tubes and strips**
Click here to search for **PCR plates**

PCR consumables

Provide your PCR reactions with the perfect housing, and take advantage of our huge choice of PCR consumables. Injection-moulded from virgin polypropylene, and quality tested thoroughly, VWR Collection PCR tubes, strips and plates come in all shapes and formats to fit your needs, as well as our adhesive sealing films and foils, available in variety of materials to suit your PCR/qPCR applications.



Compatibility Table : PCR plate with Cyclers

	Standard Profile									Low Profil									
	non-skirted		Semi Skirted							non-skirted		Semi Skirted					Full Skirted		
	Natural	White	Natural				White			Natural	White	Natural		White			Natural	White	
	VWR1732-3753	VWR1732-3754	VWR1732-3755	VWR1732-3756	VWR1732-3757	VWR1732-3758	VWR1732-3759	VWR1732-3757	VWR1732-3760	VWR1732-3761	VWR1732-3762	VWR1732-3763	VWR1732-3764	VWR1732-3765	VWR1732-3766	VWR1732-3767	VWR1732-3768	VWR1732-3769	VWR1732-3770
ABI Applied Biosystems																			
7500	•	•	•	•			•	•	•	•	•	•	•						
7700	•	•	•	•			•	•	•	•									
7500 Fast	•	•									•	•	•						
7900 HT	•	•	•	•			•	•	•	•		•	•	•					
MiniAmp			•	•			•	•	•	•									
MiniAmp Plus			•	•			•	•	•	•									
ProFlex				•	•		•	•	•	•									
Quant Studio™ 96-well			•	•			•	•	•	•									
Quant Studio™ 96-well fast										•	•	•	•	•					
SimpliAmp			•	•			•	•	•	•									
Step One Plus			•	•			•	•	•	•		•	•	•					
Veriti 0.1 ml							•	•	•	•		•	•	•					
Veriti 0.2 ml			•	•			•	•	•	•									
ViiATM 7			•	•			•	•	•	•	•	•	•	•					
Agilent																			
AriaMx										•	•						•	•	•
Mx3000P	•	•																	
Mx3005P	•	•																	
Analytik Jena / Biometra																			
qTower																	•	•	•
Tadvanced	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
TOne	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
TRobot	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
Axygen																			
MaxyGene II	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
Bio-Rad																			
C1000/C1000 Touch	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
CFX Connect			•	•			•	•	•	•	•	•	•	•			•	•	•
CFX Touch			•	•			•	•	•	•	•	•	•	•			•	•	•
S1000	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
T100	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
Eppendorf																			
Mastercycler®	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
Mastercycler® ep	•	•								•	•						•	•	•
Mastercycler® ep Realplex	•	•								•	•	•	•	•			•	•	•
Mastercycler® Gradient	•	•	•	•			•	•	•	•									
Mastercycler® Nexus	•	•	•	•			•	•	•	•									
Roche																			
LightCycler 480																	•	•	
LightCycler 96																	•	•	
Takara																			
TP 3000	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
VWR																			
UNO96	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
UNO96 HPL	•	•	•	•			•	•	•	•	•	•	•	•			•	•	•
XT96	•	•	•	•			•	•	•	•	•	•	•	•					



PCR plate 96W semi skirted 0,2 ml

PCR/QPCR PLATES, 96-WELLS, PLASTIC

PCR Plates are available in variety of styles to fit most thermal cyclers in the market.

- Broad thermal cycler compatibility
- High-contrast, black alphanumeric indexing for easy well identification
- Crystal-clear wells to easily verify sample volume
- Ultra-thin, uniform wells with optimal heat transfer and higher reaction efficiency
- Plate sealing options like film, foil, silicone mats and strip caps which reduce sample evaporation

Each plate is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.

Manufactured from high-clarity, ultra-pure polypropylene, meets FDA 21 CFR 177.1520 & USP Class VI. Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.

Bar-coded semi-skirted and full-skirted plates available.

Description	No. of wells	Profile	Skirt	Well volume (ml)	Well colour	Frame colour	Barcoded	Pk	Cat. No.
Standard profile									
PCR plates, cut corner A12	96	Standard	Non-skirted	0,2	Natural	Natural	No	50	732-3753
PCR plates, cut corner A12	96	Standard	Non-skirted	0,2	White	White	No	50	732-3754
PCR plates, cut corner A12	96	Standard	Semi-skirted	0,2	Natural	Natural	No	50	732-3755
PCR plates, cut corner A12	96	Standard	Semi-skirted	0,2	Natural	Natural	Yes	50	732-3756
PCR plates, cut corner A12	96	Standard	Semi-skirted	0,2	White	White	No	50	732-3757
PCR plates, cut corner A12	96	Standard	Semi-skirted, raised rim	0,2	Natural	Natural	No	50	732-3758
PCR plates, cut corner A12	96	Standard	Semi-skirted, raised rim	0,2	Natural	Natural	Yes	50	732-3759
PCR plates, cut corner A12	96	Standard	Semi-skirted, raised rim	0,2	White	White	No	50	732-3760
Low profile									
PCR plates, cut corner H1	96	Low	Fully skirted	0,1	Natural	Natural	No	50	732-3768
PCR plates, cut corner H1 96	96	Low	Fully skirted	0,1	Natural	Natural	Yes	50	732-3769
PCR plates, cut corner H1 96	96	Low	Fully skirted	0,1	White	White	No	50	732-3770
PCR plates, cut corner H12 96	96	Low	Non-skirted	0,1	Natural	Natural	No	50	732-3761
PCR plates, cut corner H12 96	96	Low	Non-skirted	0,1	White	White	No	50	732-3762
PCR plates, cut corner A1, ABI Fast Block	96	Low	Semi-skirted	0,1	Natural	Natural	No	50	732-3763
PCR plates, cut corner A1, ABI Fast Block	96	Low	Semi-skirted	0,1	Natural	Natural	Yes	50	732-3764
PCR plates, cut corner A1, ABI Fast Block	96	Low	Semi-skirted	0,1	White	White	No	50	732-3765
PCR plates, cut corner H12 ,Roche LC	96	Low	Semi-skirted	0,1	White	White	No	50	732-3766
PCR plates, cut corner H12, Roche LC	96	Low	Semi-skirted	0,1	White	White	Yes	50	732-3767



PCR tubes and strips

STRIPS OF PCR TUBES AND CAPS

These PCR tube strips and caps are made of PP resin and are designed to fit most popular brands of thermal cyclers.

- Certified free from DNase, RNase and human DNA
- Manufactured from polypropylene under strict quality control guidelines

Available with domed or flat caps.

Description	Colour	Capacity (ml)	Pk	Cat. No.
PCR tubes				
8-tube strips for PCR, with detached domed cap strip	Clear	0,2	250	732-1521
8-tube strips for PCR, with hinged domed cap strip	Clear	0,2	125	732-0546
8-tube strips for PCR, with individually attached domed caps, break-apart	Clear	0,2	120	732-0545
8-tube strips for PCR, with individually attached, optically clear, flat caps	Assorted	0,2	120	732-3609
8-tube strips for PCR, with separate domed cap strip	Clear	0,2	125	732-3485
12-tube strips for PCR, with separate domed 12-cap strip	Clear	0,2	80	732-0554
12-tube strips for PCR, without caps	Clear	0,2	80	732-0552
qPCR tubes				
4-tube strips, for Qiagen Rotorgene, with caps	Natural	-	250	732-1506
8-tube strips for qPCR, low profile, with individually attached, optically clear, flat caps	Clear	0,1	120	211-0339
8-tube strips for qPCR, low profile, with individually attached, optically clear, flat caps and opaque white wells	White	0,1	120	731-0367
8-tube strips for qPCR, low-profile, with individually attached, optically clear, flat caps, break-apart	Clear	0,1	120	732-3499
8-tube strips for qPCR, with individually attached, optically clear, flat caps	Clear	0,2	120	732-3608
8-tube strips for qPCR, with attached, optically clear, hinged flat cap strip	Clear	0,2	125	211-0381
8-tube strips for qPCR, with separate flat 8-cap strip	Clear	0,2	125	732-3610
8-tube strips for qPCR, without caps	White	0,2	125	732-3391
8-tube strips for qPCR/PCR, with detached flat cap strip	Clear	0,2	250	732-1520
8-tube strips for qPCR/PCR, with individually attached, optically clear, flat caps, break-apart	Clear	0,2	120	211-0338
8-tube strips for qPCR/PCR, without caps	Clear	0,2	125	732-1517
8-tube strips for qPCR/PCR, without caps	Clear	0,2	125	732-3484
Cap strips only				
8-cap strips, domed	Clear		125	732-0550
8-cap strips, domed, for 0,2 ml tube strips	Clear		125	732-1518
8-cap strips, flat, for qPCR, fitting tube strip 732-3391 & 732-1517	Clear		125	732-1519
12-cap strips, domed	Clear		80	732-0553



PCR tubes and strips

PCR TUBES, 0,2 ML

These PCR tubes with caps are made of PP resin and are designed to fit most popular brands of thermal cyclers.

- Certified free from DNase, RNase and human DNA
- Manufactured from polypropylene under strict quality control guidelines

Available with domed or flat caps.

Environmentally preferable attribute		Low Manufacturing Impact, Sustainable Packaging		
Description	Colour	Capacity (ml)	Pk	Cat. No.
Individual PCR tubes, with attached flat caps	Clear	0,2	1.000	732-0548
Individual PCR tubes, with attached domed caps	Clear	0,2	1.000	732-0547

PCR PLATES, 96-WELL, SEMI-SKIRTED

PCR Plates are available in variety of styles to fit most thermal cyclers in the market.

- Broad thermal cycler compatibility
- High-contrast, black alphanumeric indexing for easy well identification
- Crystal-clear wells to easily verify sample volume
- Ultra-thin, uniform wells with optimal heat transfer and higher reaction efficiency
- Plate sealing options like film, foil, silicone mats and strip caps which reduce sample evaporation

Each plate is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.

Manufactured from high-clarity, ultra-pure polypropylene, meets FDA 21 CFR 177.1520 & USP Class VI. Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.

Bar-coded semi-skirted and full-skirted plates available.

Description	No. of wells	Profile	Skirt	Well volume (ml)	Well colour	Frame colour	Barcoded	Pk	Cat. No.
ABI® type, corner- A12	96	Standard	Semi-skirted	0,2	Crystal clear	Clear	Yes	100	211-0283
ABI® type, corner- A12	96	Low	Semi-skirted, raised rim	0,2	Crystal clear	Clear	Yes	10	732-3487
ABI® type, corner- A12	96	Low	Semi-skirted, raised rim	0,2	Crystal clear	Assorted	No	10	732-3488
ABI® type, corner- A12	96	Low	Semi-skirted, raised rim	0,1	Crystal clear	Clear	No	10	732-3489
FAST® type, corner- A1	96	Low	Semi-skirted, raised rim	0,1	Crystal clear	Clear	No	10	732-3491
LightCycler® type	96	Low	Semi-skirted, raised rim	0,1	Crystal clear	Clear	Yes	10	732-3495
LightCycler® type	96	Low	Semi-skirted, raised rim	0,1	Crystal clear	White	Yes	10	732-3496
ABI® type, corner- A12	96	Low	Semi-skirted, raised rim	0,2	Crystal clear	Clear	Yes	10	732-3490
Straight side, corner- A12	96	Low	Semi-skirted, raised rim	0,2	Crystal clear	Clear	Yes	10	732-3494

PCR PLATES, 96-WELL, NON-SKIRTED

PCR Plates are available in variety of styles to fit most thermal cyclers in the market.

- Broad thermal cycler compatibility
- High-contrast, black alphanumeric indexing for easy well identification
- Crystal-clear wells to easily verify sample volume
- Ultra-thin, uniform wells with optimal heat transfer and higher reaction efficiency
- Plate sealing options like film, foil, silicone mats and strip caps which reduce sample evaporation

Each plate is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.

Manufactured from high-clarity, ultra-pure polypropylene, meets FDA 21 CFR 177.1520 & USP Class VI. Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.

Bar-coded semi-skirted and full-skirted plates available. Non-skirted item 3410-00 is also available in 24 and 48 well sizes.

No. of wells	Profile	Skirt	Well colour	Frame colour	Pk	Cat. No.
96	Low	Non-skirted	Clear	Clear	20	732-3486

PCR PLATES, 96-WELL

PP. These PCR plates are compatible with most thermal cyclers.

- Smooth, thin, uniform well walls ensure accurate thermal transfer
- Plates are thin and flexible
- Certified free from DNase, RNase and human genomic DNA
- Printed alphanumeric labelling and cut corner simplifies plate orientation and sample identification




732-2387



732-2388

Description	No. of wells	Profile	Skirt	Well volume (ml)	Well colour	Barcoded	Pk	Cat. No.
With coloured lettering								
PCR plate	96	Standard	Non-skirted	0,20	Clear	No	100	732-2387
PCR plate, recommended for automation	96	Standard	Semi-skirted	0,20	Clear	Yes	100	732-2390
PCR plate	96	Low	Non-skirted	0,15	Clear	No	100	732-2386
PCR plate, recommended for automation	96	Low	Semi-skirted	0,15	Clear	Yes	100	732-2388
PCR plate, recommended for automation	96	Low	Semi-skirted, raised rim	0,15	Clear	Yes	100	732-2389
PCR plate	96	Standard	Non-skirted	0,20	White	No	100	732-3388
PCR plate	96	Low	Non-skirted	0,15	White	No	100	732-3387
PCR plate, recommended for automation	96	Low	Semi-skirted	0,15	White	Yes	100	732-3389
PCR plate, recommended for automation	96	Low	Semi-skirted, raised rim	0,15	White	Yes	100	732-3390
PCR plate	96	Low	Fully skirted	0,10	Purple	No	100	211-0302
Without coloured lettering								
PCR plate	96	Standard	Non-skirted	0,20	Clear	No	100	211-0262
PCR plate, raised well	96	Standard	Non-skirted	0,20	Clear	No	100	211-0269
PCR plate	96	Low	Fully skirted	0,10	Clear	No	100	211-0297
PCR plate	96	Standard	Non-skirted	0,20	Assorted	No	100	211-0263
PCR plate	96	Standard	Non-skirted	0,20	Blue	No	100	211-0264
PCR plate	96	Standard	Non-skirted	0,20	Green	No	100	211-0265
PCR plate	96	Standard	Non-skirted	0,20	Purple	No	100	211-0266
PCR plate	96	Standard	Non-skirted	0,20	Red	No	100	211-0267
PCR plate	96	Standard	Non-skirted	0,20	Yellow	No	100	211-0268
PCR plate	96	Standard	Fully skirted	0,20	Assorted	No	100	211-0298
PCR plate	96	Low	Fully skirted	0,10	Blue	No	100	211-0300
PCR plate	96	Standard	Fully skirted	0,20	Green	No	100	211-0301
PCR plate	96	Low	Fully skirted	0,10	Red	No	100	211-0303
PCR plate	96	Standard	Fully skirted	0,20	Yellow	No	100	211-0304

Description	Colour	Pk	Cat. No.
Cap strips only			
8-cap strips, domed	Clear	125	732-0550
8-cap strips, domed, for 0,2 ml tube strips	Clear	125	732-1518
8-cap strips, flat, for qPCR, fitting tube strip 732-3391 & 732-1517	Clear	125	732-1519



qPCR PLATES, 96-WELL

PP. These white qPCR plates and optically-clear cap closures are suitable for Real-Time qPCR applications. White qPCR plates are designed to enable sensitive and accurate fluorescence detection. When used together with ultra-clear caps or optical seals, these products will increase sensitivity and reduce variability in qPCR assays.

- Smooth, thin, uniform well walls ensure accurate thermal transfer
- Wells are slightly raised to accept optically-clear strip caps or sealing film
- Certified free from DNase, RNase and human genomic DNA

Working capacity: 200 µl

No. of wells	Profile	Skirt	Well volume (ml)	Well colour	Frame colour	Pk	Cat. No.
96	Standard	Non-skirted	0,2	White	White	100	211-0313
96	Low	Fully skirted	0,2	White	White	100	211-0315
96	Standard	Semi-skirted	0,2	White	White	100	211-0317

Description	Pk	Cat. No.
Accessories		
Optically clear 8-cap strips for Real-Time PCR plates	125	211-0350

PCR PLATES, 384-WELL, SKIRTED

PCR Plates are available in variety of styles to fit most thermal cyclers in the market.

- Broad thermal cycler compatibility
- High-contrast, black alphanumeric indexing for easy well identification
- Crystal-clear wells to easily verify sample volume
- Ultra-thin, uniform wells with optimal heat transfer and higher reaction efficiency
- Plate sealing options like film, foil, silicone mats and strip caps which reduce sample evaporation

Each plate is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.

Manufactured from high-clarity, ultra-pure polypropylene meets FDA 21 CFR 177.1520 & USP Class VI. Certified free of detectable DNase, RNase, DNA, PCR inhibitors, and tested pyrogen-free.

Bar-coded semi-skirted and full-skirted plates available.

Environmentally preferable attribute	Low Manufacturing Impact, Sustainable Packaging
---	---

Description	No. of wells	Skirt	Well colour	Frame colour	Barcoded	Pk	Cat. No.
ABI® type, one-notch, corner- A24	384	Fully skirted	Crystal clear	Clear	No	10	732-3422



PCR PLATES, 384-WELL

PP. These PCR plates are compatible with most thermal cyclers, and are ideal for high throughput screening thermal cycler applications.

- Smooth, thin, uniform well walls ensure accurate thermal transfer
- Wells are slightly raised to accommodate sealing mats, films or foils
- Plates are skirted to allow barcoding, and include a frosted labelling area
- Lot tested and certified free from DNase, RNase and human genomic DNA
- Printed alphanumeric labelling simplifies plate and sample identification

Working capacity: 25 µl

No. of wells	Skirt	Well colour	Frame colour	Pk	Cat. No.
384	Fully skirted	Clear	Clear	100	211-0305

Selection Guide by Product Application : Sealing Films

Product Description	EU Catalog No.	Non-Skirted / Raised-Rim Microplates	Flat-Top Microplates	Dual Component (Polycarbonate/ Polypropylene)	PCR 8-Tube Strips / Strip Wells/ Single Plate Row	PCR	qPCR (Real-Time)	Short-Term Storage/Incubation/ELISA	Long-Term Storage	Cell/Tissue/Bacterial Culture	Protein Crystallography	HPLC	Fluorescence/Luminescence	Sonication
VWR® Heat-Resistant Films for Real-Time qPCR, Ultra-Clear Polyester, Raised Rim	391-1295	•	•	•			•							
VWR® Heat-Resistant Films for Real-Time qPCR, Ultra-Clear Polyester	391-1258		•	•			•							
VWR® Polyolefin Films with Silicone Adhesive for qPCR, Storage, and Protein Crystallization	391-0189	•	•				•		•		•	•		
VWR® Polyolefin Silicone Film	391-0624	•	•				•		•		•	•		
VWR® Heat-Resistant Strip Films for PCR	731-0321		•			•								
VWR® Heat-Resistant Polypropylene Film for Raised-Rim Plates	391-1294	•				•								
VWR® Ultimate Seal Films for PCR and Microplates	391-0622	•	•			•								
VWR® Heat-Resistant Films for PCR, Polypropylene	391-1254		•			•								
VWR® Heat-Resistant Films for PCR, Polypropylene	391-1255		•			•								
VWR® Aluminum Foils for PCR and Cold Storage	391-1256		•	•		•			•					
VWR® Aluminum Foils for PCR and Cold Storage	391-1257		•	•		•			•					
VWR® Aluminum Foils for 96- and 384-Well Plates	391-1281		•			•			•					
VWR® Aluminum Foils for 96- and 384-Well Plates	391-1282	•				•			•					
VWR® Precut Pierceable Vinyl Films for Robotics	391-1286		•									•		•
VWR® Precut Pierceable Vinyl Films for Robotics	391-1287		•									•		•
VWR® Clear-Zone Pierceable Films for Robotics	391-1264		•									•		
VWR® Adhesive Sealing Films in Roll Format for Automated Sealers	731-0307		•			•			•					
VWR® Adhesive Sealing Films in Roll Format for Automated Sealers	731-0308		•				•		•		•	•		
Adhesive sealing film qPCR/storage	391-0644	•	•				•		•		•	•		
Adhesive sealing foil strips pierc. PCR	391-0637		•		•	•								
Adhesive sealing film strips PE PCR	391-0638		•		•	•								
Adhesive sealing foil roll piercab. PCR	391-0641	•	•			•								
Adhesive sealing film PCR/storage large	391-0643		•			•								




Solutions for discovery

VWR® electrophoresis

Instruments, reagents and accessories for

- Nucleic acid electrophoresis
- Protein electrophoresis
- Blotting
- Bio imaging



For your copy visit vwr.com

Selection Guide by Product Attribute : Sealing Films

Product Description	EU Catalog No.	Functional Temperature Range	Sterile	Rnase/Dnase Free	Chemically Resistant	Pierceable	Optically Clear	Transparent	Light Blocking	Breathable	Film / Foil Thickness	Adhesive, Thickness
VWR® Heat-Resistant Films for Real-Time qPCR, Ultra-Clear Polyester, Raised Rim	391-1295	-40°C to +120°C		●			●				Polyester 50.8µm 2mil	Acrylic 29.2µm 1.15mil
VWR® Heat-Resistant Films for Real-Time qPCR, Ultra-Clear Polyester	391-1258	-40°C to +120°C		●			●				Polyester 50.8µm 2mil	Acrylic 29.2µm 1.15mil
VWR® Polyolefin Films with Silicone Adhesive for qPCR, Storage, and Protein Crystallization	391-0189	-80°C to +105°C		●	●		●				Polyolefin 50.8µm 2mil	Silicone 50.8µm 2mil
VWR® Polyolefin Silicone Film	391-0624	-80°C to +105°C	●	●	●		●				Polyolefin 50.8µm 2mil	Silicone 50.8µm 2mil
VWR® Heat-Resistant Strip Films for PCR	731-0321	-40°C to +120°C		●				●			Polypropylene 50.8µm 2mil	Acrylic 26.7µm 1.05mil
VWR® Heat-Resistant Polypropylene Film for Raised-Rim Plates	391-1294	-40°C to +120°C		●				●			Polyester 50.8µm 2mil	Acrylic 41.9µm 1.65mil
VWR® Ultimate Seal Films for PCR and Microplates	391-0622	-20°C to +120°C		●				●			Polyester 50.8µm 2mil	Acrylic 41.9µm 1.65mil
VWR® Heat-Resistant Films for PCR, Polypropylene	391-1254	-40°C to +120°C		●				●			Polypropylene 50.8µm 2mil	Acrylic 24.1µm .95mil
VWR® Heat-Resistant Films for PCR, Polypropylene	391-1255	-40°C to +120°C		●				●			Polypropylene 50.8µm 2mil	Acrylic 24.1µm .95mil
VWR® Aluminum Foils for PCR and Cold Storage	391-1256	-80°C to +120°C		●	●	●			●		Aluminum foil 35.6µm 1.4mil	Acrylic 27.9µm 1.1mil
VWR® Aluminum Foils for PCR and Cold Storage	391-1257	-80°C to +120°C		●	●	●			●		Aluminum foil 35.6µm 1.4mil	Acrylic 27.9µm 1.1mil
VWR® Aluminum Foils for 96- and 384-Well Plates	391-1281	-80°C to +120°C		●		●			●		Aluminum foil 38.1µm 1.5mil	Acrylic 37.6µm 1.48mil
VWR® Aluminum Foils for 96- and 384-Well Plates	391-1282	-80°C to +120°C		●		●			●		Aluminum foil 38.1µm 1.5mil	Acrylic 37.6µm 1.48mil
VWR® Precut Pierceable Vinyl Films for Robotics	391-1286	-40°C to +90°C					●	●			Vinyl 88.9µm 3.5mil	Acrylic 17.8µm 0.7mil
VWR® Precut Pierceable Vinyl Films for Robotics	391-1287	-40°C to +90°C	●				●	●			Vinyl 88.9µm 3.5mil	Acrylic 17.8µm 0.7mil
VWR® Clear-Zone Pierceable Films for Robotics	391-1264	-40°C to +90°C					●	●			Polyethylene 69.9µm 2.75mil	Acrylic 26.7µm 1.05mil
VWR® Adhesive Sealing Films in Roll Format for Automated Sealers	731-0307	-80°C to +120°C		●		●			●		Aluminum foil 38.1µm 1.5mil	Acrylic 37.6µm 1.48mil
VWR® Adhesive Sealing Films in Roll Format for Automated Sealers	731-0308	-80°C to +105°C		●	●		●				Polyolefin 50.8µm 2mil	Silicone 50.8µm 2mil
Adhesive sealing film qPCR/storage	391-0644	-80°C to +105°C		●	●		●				Polyolefin 50.8µm 2mil	Silicone 50.8µm 2mil
Adhesive sealing foil strips pierc. PCR	391-0637	-20°C to +120°C		●		●			●		Aluminum Foil 50.8µm 2mil	Acrylic 36.8µm 1.45mil
Adhesive sealing film strips PE PCR	391-0638	-20°C to +120°C		●				●			Polypropylene 50.8µm 2mil	Acrylic 41.9µm 1.65mil
Adhesive sealing foil roll piercab. PCR	391-0641	-20°C to +120°C		●		●			●		Aluminum Foil 50.8µm 2mil	Acrylic 36.8µm 1.45mil
Adhesive sealing film PCR/storage large	391-0643	-20°C to +120°C		●				●			Polyester 50.8µm 2mil	Acrylic 41.9µm 1.65mil



ADHESIVE PCR FILM SEALS

These heat resistant 74,9 µm thick films are designed for thermal cycling applications. Polypropylene films are not pierceable. For PCR applications, where piercing with pipette tips or robotic probes is required for product recovery, use aluminium foil films. For Real-Time PCR applications, where maximum optical clarity is required, use optically-clear polyester films.


- Recommended for temperatures from -40 to +120 °C
- Certified free from DNase, RNase and nuclease

Each film LxD: 135,1x79,4 mm with sufficient sealing area for all PCR plates.

Length with end tabs removed: 123,1 mm

* stronger, thicker adhesive and cut to fit raised-rim plates

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Adhesive PCR film seals	Transparent	-	74,9	100	391-1254
Adhesive PCR film seals	Transparent	+	74,9	100	391-1255
Adhesive PCR film strips to seal 2x8 wells	Transparent	-	74,9	200	731-0321
Adhesive PCR film seals*	Transparent	-	92,7	100	391-1294



Adhesive sealing films for PCR and storage

ADHESIVE SEALING FILMS FOR PCR AND STORAGE

Polyester seals feature a strong adhesive to minimize evaporation during PCR and storage applications.

- Ideal for PCR & storage
- Certified DNase- and RNase- free
- Temperature range is -40 to +120 °C
- Non-pierceable
- Minimizes evaporation

Removes easily without leaving adhesive residue.


Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Sealing film	Transparent	-	92,7	100	391-0622
Sealing film, large	Transparent	-	92,7	100	391-0643

ALUMINUM SEALING STRIPS FOR 8-STRIP PCR TUBES

Easier to apply and remove as compared to standard molded PCR strip caps.

- Ideal for PCR
- User-friendly alternative to strip-caps; no sore thumbs and fingers guaranteed
- Seals a single 8-well row of a 96-well plate or a single PCR 8-tube strip
- Certified DNase- and RNase-free
- 2.0 mil pierceable foil
- Recommended temperature range from -20 to +120 °C

6 strips are precut on a continuous liner for ease of use.



Aluminum sealing strips for 8-strip PCR tubes

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
1x8 sealing strips	Light blocking	-	87,6	300	391-0637

ADHESIVE SEALING STRIPS FOR 8-STRIP PCR TUBES

Easier to apply and remove as compared to standard molded PCR strip caps.




- Ideal for PCR and storage
- Alternative to strip-caps; no sore thumbs and fingers guaranteed
- Seals a single 8-well row of a 96-well plate or a single PCR 8-tube strip
- Certified DNase- and RNase-free
- 2.0 mil transparent PE film
- Recommended temperature range from -20 to +120 °C

8 strips are precut on a continuous liner for ease of use.

Adhesive sealing strips for 8-strip PCR tubes

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
1x8 sealing strips	Transparent	-	92,7	400	391-0638



ADHESIVE PCR FOIL SEALS

These soft, non permeable 63,5 µm thick aluminium foils, with strong medical grade adhesive, eliminate the need for heat sealing devices or mats during thermal cycling. Compared to other aluminium foils, these foils have less tendency to roll back on themselves when removing the backing paper, and fit well to the plate during application.

- Recommended for temperatures from -80 to +120 °C
- Easily pierceable with pipette tips and robotic probes
- Excellent vapour barrier, virtually no sample evaporation
- Certified free from DNase, RNase and nuclease

Each foil LxD: 142,9x82,6 mm with sufficient sealing area for all PCR plates.

Length with end tabs removed: 125,4 mm

Sterile product is packed in tamper-evident bags of 25/bag.

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Aluminium foils	Light blocking	+	63,5	50	391-1257
Aluminium foils	Light blocking	-	63,5	100	391-1256




ALUMINIUM FOIL SEALS FOR PCR AND STORAGE (96- AND 384-WELL PLATES)

Aluminium foils, 75,7 µm thick, for use with 96- or 384-well plates. Fit inside the rim of raised rim plates. These foils have one partial-width end tab. Available non sterile only.

- Recommended for temperatures from -80 to +120 °C
- Certified free from DNase, RNase and nucleic acids

Each foil LxD: 127,0x77,8 mm, including single 9,5 mm end tab.

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Aluminium foils for 96-well plates	Light blocking	-	75,7	100	391-1282
Aluminium foils for 384-well plates	Light blocking	-	75,7	100	391-1281



Adhesive sealing foil roll

ALUMINIUM FOIL SEALS FOR PCR

Precut aluminum foil seals on a roll for PCR applications. Supplied on a roll with holder for ease of use.

- Ideal for PCR
- Rapidly seal microplates; faster than applying single sheeted film
- 2 mil pierceable foil
- Certified DNase- and RNase-free
- Recommended temperature range from -20 to +120 °C

Self-standing, non-slip dispenser is supplied for ease of use.

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Sealing foil, roll format	Light blocking	-	87,6	500	391-0641



SEALING FILM FOR AUTOMATION

Adhesive sealing film or foil rolls constructed on plastic cores, for use with high throughput automated microplate sealers.

- Continuous rolls - no splices
- Film or foil does not extend beyond edge of plastic core
- Robust static-free packaging provides protection during shipment and facilitates accurate alignment and adhesion

Sealing film for qPCR, sitting-drop protein crystallisation, HTS and cold storage

Clear polyolefin with ultra-strong, pressure-sensitive silicone adhesive for qPCR, sitting-drop protein crystallisation, HTS, storage and general incubation applications.


- DMSO resistant
- Nuclease-, nucleic acid- and pyrogen-free
- Provides as strong a bond to the plate as heat sealing films without the drawbacks of heat

Sealing foil for PCR, HTS and cold storage

Pierceable aluminium foil with pressure-sensitive acrylic adhesive for PCR, HTS and cold storage.

- DMSO resistant
- Easily pierceable for sample recovery
- Nuclease-, nucleic acid- and pyrogen-free

Description	Optical property	Sterile	Temp. range (°C)	Thickness (µm)	Pk	Cat. No.
Sealing film, polyolefin, DNase- and RNase-free	Optically clear	-	-80...+105	101,6	1 Roll	731-0308
Sealing foil, aluminium	Light blocking	-	-80...+120	75,7	1 Roll	731-0307



ULTRA-CLEAR FILMS FOR QPCR


Transparent polyester films with strong, non absorbing, non fluorescing, medical grade adhesive for superior performance in qPCR applications. Supplied non sterile.

- Recommended for temperatures from -40 to +120 °C
- Ultra-high optical clarity
- Certified free from DNase, RNase and nuclease

Each film LxD: 142,9x79,4 mm

Length with end tabs removed: 121,9 mm

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Optically-clear	Optically clear	-	80,0	100	391-1258
Optically-clear, for raised rim plates	Optically clear	-	80,0	100	391-1295




SEALING FILMS, NON-TACKY, FOR QPCR, STORAGE AND CRYSTALLISATION

50 µm polyolefin films with 50 µm inert encapsulated silicone adhesive. Especially suited for qPCR, storage and protein crystallisation applications. The encapsulated silicone adhesive is non tacky until pressed against the sealing surface, at which time adhesive is released only in sealing areas to form the strongest available heat resistant seal around each well on the plate.

- High optical clarity, minimal to no autofluorescence
- Chemically inert; no extractables except at extreme pH; DMSO resistant for HTS
- Heat resistant, recommended for temperatures from -80 to +105 °C
- Certified free from DNase, RNase, and nucleic acid

Sized to fit within the edges of raised-rim 96-well plates (76,2x133,4 mm). Two end tabs assist in positioning the film on the plate.

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Sealing film	Optically clear	-	101,6	100	391-0189
Sealing film	Optically clear	+	101,6	50	391-0624
Sealing film	Optically clear	-	101,6	25	391-0644




PIERCEABLE FILMS FOR ROBOTICS, CLEAR-ZONE

With a polyethylene top layer and inert white polypropylene and acrylic adhesive sublayer, this easy-piercing sealing film allows direct sample recovery with single or multichannel pipettors and robotic probes.

- Clear zone above each well is free of adhesive
- Minimise sample exposure to adhesive
- Recommended temperature range -40 to +90 °C

The adhesive-free zone above each well prevents adhesive fouling of the pipette tip or robotic probe, minimising the amount of adhesive that comes into contact with the sample. Suitable for 96-well plates only.

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Clear-Zone films	Transparent	-	96,6	50	391-1264



PRECUT PIERCEABLE VINYL FILMS FOR ROBOTICS

These 106,7 µm thick vinyl films with 20 µm adhesive layer are designed for temporary protection of samples in 96-well plates from contamination and evaporation.

- Protect samples and limit evaporation short-term
- Precut flaps bend inward without fouling probes or tips
- Flaps close for continued protection after sampling
- Recommended temperature range -40 to +90 °C

A precut pattern over each well separates the film into four flaps that bend inward easily when pushed by a robotic probe or pipette tip, allowing access to the sample without coring or adhesive fouling. The resilient flaps regain their original position after sampling for continued sample protection. For long-term sample protection after sampling, a continuous film should be applied as a second layer. Suitable for sealing all standard 96-well plates.

Each film LxD: 145,5x79,4 mm

Description	Optical property	Sterile	Thickness (µm)	Pk	Cat. No.
Precut vinyl films	Transparent	-	106,7	100	391-1286
Precut vinyl films	Transparent	+	106,7	50	391-1287



PCR RACKS, REVERSIBLE

- PCR side of the rack has 168 wells that can hold 8-well or 12-well tube strips, or one 0,2 ml tube per well
- Opposite side of the rack has 40 wells that can hold 0,5 ml tubes
- Both sides of the rack have 12 wells that can hold 1,5 ml tubes
- Simple to use, removable hinged lid snaps in place on either side of the rack

Designed to hold 0,2; 0,5 or 1,5 ml tubes.

Assorted pack includes one each of blue, green, purple, yellow and orange.

Description	Colour	Pk	Cat. No.
PCR racks	Assorted colours	5	211-0198



Microcentrifuge tubes

MICROCENTRIFUGE TUBES, HIGH PERFORMANCE

These tubes are made from USP Class VI PP, which ensures a clear view of the pellet. They incorporate a high strength wall, and can be centrifuged to twice the g-force that most other microcentrifuge tubes can withstand.

- Amber coloured for light-sensitive substances
- Ideal for vortexing and shaking tables
- RNase-, DNase- and endotoxin-free
- Engraved graduations ensure accuracy
- Fit all common rotors
- Ideal for phenol/chloroform extractions

These tubes also feature flat, pierceable, frosted caps and a side labelling surface. The ergonomically designed caps, with smooth surfaces, help reduce the potential for repetitive strain injuries.

Steam autoclavable at 121 °C and freezable to -80 °C.

Manufactured in a Class 100000 cleanroom environment

Manufactured under ISO 13485:2016 and ISO 9001:2015 quality management systems

Assorted pack contains 100 tubes each of green, blue, red, orange and yellow in separate bags.

Because tubes have a tightly fitting cap, they should not be used for boiling.

Capacity (ml)	Type	Colour	Version	RCF max. (xg)	Pk	Cat. No.
0,5 ml tubes						
0,5	With attached lid	Natural	Graduated, conical bottom, non sterile	35000	1.000	525-1172
0,5	With attached lid	Assorted	Graduated, conical bottom, non sterile	35000	1.000	525-1170
0,5	With attached lid	Blue	Graduated, conical bottom, non sterile	35000	1.000	525-1171
0,5	With attached lid	Green	Graduated, conical bottom, non sterile	35000	1.000	525-1177
0,5	With attached lid	Orange	Graduated, conical bottom, non sterile	35000	1.000	525-1173
0,5	With attached lid	Purple	Graduated, conical bottom, non sterile	35000	1.000	525-1174
0,5	With attached lid	Red	Graduated, conical bottom, non sterile	35000	1.000	525-1175
0,5	With attached lid	Yellow	Graduated, conical bottom, non sterile	35000	1.000	525-1176
1,5 ml tubes						
1,5	With attached lid	Natural	Graduated, conical bottom, non sterile	40000	500	525-1164
1,5	With attached lid	Amber	Graduated, conical bottom, non sterile	40000	500	525-1161
1,5	With attached lid	Assorted	Graduated, conical bottom, non sterile	40000	500	525-1167
1,5	With attached lid	Blue	Graduated, conical bottom, non sterile	40000	500	525-1162
1,5	With attached lid	Green	Graduated, conical bottom, non sterile	40000	500	525-1163
1,5	With attached lid	Orange	Graduated, conical bottom, non sterile	40000	500	525-1165
1,5	With attached lid	Purple	Graduated, conical bottom, non sterile	40000	500	525-1166
1,5	With attached lid	Red	Graduated, conical bottom, non sterile	40000	500	525-1168
1,5	With attached lid	Yellow	Graduated, conical bottom, non sterile	40000	500	525-1169



Microcentrifuge tubes

DISPOSABLE MICROCENTRIFUGE TUBES

These tubes are made from USP Class VI PP.

- Amber coloured for light-sensitive substances
- RNase-, DNase- and endotoxin-free
- Engraved graduations ensure accuracy
- Caps will not pop open during centrifuging
- Caps can be penetrated by a syringe

These tubes also feature flat, pierceable, frosted caps and a side labelling surface.

Autoclavable at 121 °C and freezable to -80 °C.

Manufactured in a Class 100000 cleanroom environment

Manufactured under ISO 13485:2016 and IS O 9001:2015 quality management systems

Assorted pack contains 100 tubes each of red, orange, yellow, green and blue.

Capacity (ml)	Type	Colour	Version	RCF max. (×g)	Pk	Cat. No.
0,5 ml tubes						
0,5	With attached cap	Natural	Graduated, conical bottom, non sterile	20000	1.000	525-1220
0,5	With attached cap	Amber	Graduated, conical bottom, non sterile	20000	1.000	525-1246
0,5	With attached cap	Assorted	Graduated, conical bottom, non sterile	20000	1.000	525-1222
0,5	With attached cap	Blue	Graduated, conical bottom, non sterile	20000	1.000	525-1250
0,5	With attached cap	Green	Graduated, conical bottom, non sterile	20000	1.000	525-1251
0,5	With attached cap	Orange	Graduated, conical bottom, non sterile	20000	1.000	525-1248
0,5	With attached cap	Purple	Graduated, conical bottom, non sterile	20000	1.000	525-1249
0,5	With attached cap	Red	Graduated, conical bottom, non sterile	20000	1.000	525-1221
0,5	With attached cap	Yellow	Graduated, conical bottom, non sterile	20000	1.000	525-1247
1,5 ml tubes						
1,5	With attached cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1126
1,5	With attached cap	Natural	Graduated, conical bottom, sterile	20000	500	525-1178
1,5	With attached cap	Amber	Graduated, conical bottom, non sterile	20000	500	525-1223
1,5	With attached cap	Assorted	Graduated, conical bottom, non sterile	20000	500	525-1230
1,5	With attached cap	Blue	Graduated, conical bottom, non sterile	20000	500	525-1228
1,5	With attached cap	Green	Graduated, conical bottom, non sterile	20000	500	525-1229
1,5	With attached cap	Orange	Graduated, conical bottom, non sterile	20000	500	525-1225
1,5	With attached cap	Purple	Graduated, conical bottom, non sterile	20000	500	525-1227
1,5	With attached cap	Red	Graduated, conical bottom, non sterile	20000	500	525-1226
1,5	With attached cap	Yellow	Graduated, conical bottom, non sterile	20000	500	525-1224
2,0 ml tubes						
2,0	With attached cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1136
2,0	With attached cap	Amber	Graduated, conical bottom, non sterile	20000	500	525-1127
2,0	With attached cap	Assorted	Graduated, conical bottom, non sterile	20000	500	525-1135
2,0	With attached cap	Blue	Graduated, conical bottom, non sterile	20000	500	525-1133
2,0	With attached cap	Green	Graduated, conical bottom, non sterile	20000	500	525-1134
2,0	With attached cap	Orange	Graduated, conical bottom, non sterile	20000	500	525-1129
2,0	With attached cap	Purple	Graduated, conical bottom, non sterile	20000	500	525-1131
2,0	With attached cap	Red	Graduated, conical bottom, non sterile	20000	500	525-1130
2,0	With attached cap	Yellow	Graduated, conical bottom, non sterile	20000	500	525-1128



Microcentrifuge Tubes with Socket Screw-Cap

MICROCENTRIFUGE TUBES WITH SOCKET SCREW CAPS

These clear, high grade PP tubes are ideal for sample processing applications.

- RNase-, DNase- and endotoxin-free

Excellent for freezer storage or secure sample handling. Screw caps incorporate silicone O-rings for a leakproof seal and accept colour coded inserts. Inserts snap into the top of the cap and create a flush surface. Tubes are available with standard or loop caps, which tether the cap to the tube.

Caps and tubes are steam autoclavable at 121 °C and freezable to -80 °C.

Manufactured in a Class 100000 cleanroom environment

Manufactured under ISO 13485:2016 and ISO 9001:2015 quality management systems

Sterile tubes are sterilised by gamma irradiation SAL 10⁻⁶ (ISO 11137).

Capacity (ml)	Type	Colour	Version	RCF max. (xg)	Pk	Cat. No.
Tubes with green colour caps in separate resealable bags						
0,5	With green screw cap	Natural	Non graduated, conical bottom, non sterile	20000	500	525-1231
0,5	With green screw cap	Natural	Non graduated, freestanding, non sterile	20000	500	525-1233
1,5	With green screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1236
1,5	With green screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1237
2,0	With green screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1240
2,0	With green screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1241
Sterile tubes with attached natural colour caps						
0,5	With natural screw cap	Natural	Non graduated, conical bottom, sterile	20000	500	525-0642
0,5	With natural screw cap	Natural	Non graduated, freestanding, sterile	20000	500	525-0644
1,5	With natural screw cap	Natural	Graduated, conical bottom, sterile	20000	500	525-0646
1,5	With natural screw cap	Natural	Graduated, freestanding, sterile	20000	500	525-0648
2,0	With natural screw cap	Natural	Graduated, conical bottom, sterile	20000	500	525-0650
2,0	With natural screw cap	Natural	Graduated, freestanding, sterile	20000	500	525-0652
Tubes with green colour loop style caps						
0,5	With green loop style screw cap	Natural	Non graduated, conical bottom, non sterile	20000	500	525-1234
0,5	With green loop style screw cap	Natural	Non graduated, freestanding, non sterile	20000	500	525-1235
1,5	With green loop style screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1238
1,5	With green loop style screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1239
2,0	With green loop style screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1243
2,0	With green loop style screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1244



Microcentrifuge tubes

MICROCENTRIFUGE TUBES WITH FLAT SCREW CAPS

These tubes are made from ultra-clear USP Class VI PP that allows easy viewing of contents.

- RNase-, DNase- and endotoxin-free
- Engraved graduations ensure accuracy
- Fit most common rotors
- Smooth surface for applying labels

These tubes feature universal screw threads for use with most popular brand screw caps. Low profile caps include an ethylene propylene O-ring, ensuring a leakproof seal.

Tubes and caps are autoclavable at 121 °C and freezable to -80 °C.

Manufactured in a Class 100000 cleanroom environment

Manufactured under ISO 13485:2016 and ISO 9001:2015 quality management systems

Sterile tubes are radiation sterilised SAL 10⁻⁶ (ISO 11137). Certificate available for every lot.

Tubes are packed 500 per bag. Assorted pack of caps (525-1187) contains equal quantities of red, orange, yellow, blue and green.

Capacity (ml)	Type	Colour	Version	RCF max. (xg)	Pk	Cat. No.
Non sterile microcentrifuge tubes with natural colour caps in separate resealable bags						
0,5	With natural screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1144
0,5	With natural screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1157
1,5	With natural screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1145
1,5	With natural screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1158
2,0	With natural screw cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1146
2,0	With natural screw cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1160
Sterile microcentrifuge tubes with natural colour caps attached						
0,5	With natural screw cap	Natural	Graduated, conical bottom, sterile	20000	500	525-1140
0,5	With natural screw cap	Natural	Graduated, freestanding, sterile	20000	500	525-1155
1,5	With natural screw cap	Natural	Graduated, conical bottom, sterile	20000	500	525-1141
1,5	With natural screw cap	Natural	Graduated, freestanding, sterile	20000	500	525-1159
2,0	With natural screw cap	Natural	Graduated, conical bottom, sterile	20000	500	525-1143
2,0	With natural screw cap	Natural	Graduated, freestanding, sterile	20000	500	525-1156
Microcentrifuge tubes only						
0,5	Without cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1137
0,5	Without cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1147
1,5	Without cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1138
1,5	Without cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1148
2,0	Without cap	Natural	Graduated, conical bottom, non sterile	20000	500	525-1139
2,0	Without cap	Natural	Graduated, freestanding, non sterile	20000	500	525-1149

Colour	Pk	Cat. No.
Flat screw caps for microcentrifuge tubes		
Assorted	500	525-1187
Blue	500	525-1180
Green	500	525-1179
Natural	500	525-1186
Orange	500	525-1183
Purple	500	525-1181
Red	500	525-1182
White	500	525-1185
Yellow	500	525-1184

J.T.BAKER® ROBOTIC TIPS FOR HAMILTON INSTRUMENTS

Available in conductive version, with or without filter.

- Manufactured under cleanroom conditions (cleanroom class 8)
- State of the art quality control
- Ideal for high-throughput analysis
- Barcoded trays for automated identification
- Packed in blister

Type `Hamilton` tips are compatible to the following platforms:

Hamilton Microlab Star, Nimbus, Vantage, Microlab Prep.



J.T.Baker® Conductive Tips 300 µl for Hamilton instruments - in Tray

Volume (µl)	Filter	Tip colour	Packed	Pk	Cat. No.
50		Black	12× 5 racks, 96 each, blistered	5.760	49010-0003
50	+	Black	12× 5 racks, 96 each, blistered	5.760	49010-0105
300		Black	12× 5 racks, 96 each, blistered	5.760	49008-0003
300	+	Black	12× 5 racks, 96 each, blistered	5.760	49008-0105
1000		Black	8× 5 racks, 96 each, blistered	3.840	49009-0005
1000	+	Black	8× 5 racks, 96 each, blistered	3.840	49009-0105

J.T.BAKER® ROBOTIC TIPS FOR TECAN INSTRUMENTS

Available in conductive and non-conductive version, with or without filter.

- Manufactured under cleanroom conditions (cleanroom class 8)
- State of the art quality control
- Ideal for high-throughput analysis
- Packed in cardboard or blister

Type `Tecan` is compatible with the following robotic instruments:

Tecan: Genesis, Freedom EVO, MiniPrep, Cavro, Fluent

Qiagen/Corbett: CAS 1200/4200/4800, QIAgility (49002-0106 only)

Perkin Elmer (Packard): Multiprobe II, Janus (with multiprobe head)

Abbott Diagnostic: m2000

Adaltis: NEXgen

AXA Diagnostics: SkyLAB 752

Hologic: Panther Fusion

NGNY: AQUA System



J.T.Baker® Conductive Tips 1000 µl for Tecan instruments in S-Tray

Volume (µl)	Filter	Tip colour	Packed	Pk	Cat. No.
Conductive tips					
50		Black	18× 10 racks, 96 each, cardboard, S-Tray	17.280	49006-0000
50		Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49006-0006
50		Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49006-0007
50	+	Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49006-0106
50	+	Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49006-0107
200		Black	18× 10 racks, 96 each, cardboard, S-Tray	17.280	49002-0000
200		Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49002-0006
200		Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49002-0007
200	+	Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49002-0106
200	+	Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49002-0107
1000		Black	10× 10 racks, 96 each, cardboard, S-Tray	9.600	49003-0000
1000		Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49003-0006
1000		Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49003-0007
1000	+	Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49003-0106
1000	+	Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49003-0107
Conductive tips, wide bore					
1000		Black	10× 10 racks, 96 each, cardboard, S-Tray	9.600	49003-0000
1000		Black	12× 2 racks, 96 each, blistered, S-Tray	2.304	49003-0006
1000		Black	12× 2 racks, 96 each, blistered, C-Tray	2.304	49003-0106
Clear tips					
200	+	Clear	18× 10 racks, 96 each, cardboard, S-Tray	17.280	49102-0000
1000		Clear	10× 10 racks, 96 each, cardboard, S-Tray	9.600	49103-0000



J.T.Baker® Conductive Tips 1000 µl for Roche MagNA Pure in blister packaging

J.T.BAKER® ROBOTIC TIPS FOR ROCHE MAGNA PURE

Available in conductive version, with filter, non-sterile.

- Manufactured under cleanroom conditions (cleanroom class 8)
- State of the art quality control
- Barcoded trays for automated identification
- Ideal for high-throughput analysis
- Packed in blister

The tips are designed for the use with Roche MagNA Pure instrument.

Volume (µ l)	Filter	Tip colour	Packed	Pk	Cat. No.
1000	+	Black	8× 5 racks, 96 each, blistered	3.840	49209-0105

J.T.BAKER® ROBOTIC TIPS FOR STRATEC INSTRUMENTS

Available in conductive version, with or without filter.

- Manufactured under cleanroom conditions (cleanroom class 8)
- Conductive tips for Liquid Level Sensing
- Ideal for dosing even small volumes with high precision and without contamination

Type Stratec is compatible with the following robotic instruments:

Abbott Diagnostic: en4lisa; bioMérieux: DaVinci; Bio-Rad: Evolis; Biotest: QuickStep; Dade Behring (Siemens): BEP 2000 Advance; DiaSorin: ETI-Max 3000; Euroimmun: Euroimmun Analyzer; Grifols: Triturus; Human: Elisys Quattro; QIAgen: Biorobot 3000, 8000, 9600, 9604, Biorobot universal; Robimes: Elisia Plus; Stratec: Gemini; Virion\Serion: Immunomat TWINsystem



J.T.Baker® Conductive Tips 300 µl for Stratec instruments - in Tray

Volume (µ l)	Filter	Tip colour	Packed	Pk	Cat. No.
300		Black	18× 10 racks, 96 each, cardboard	17.280	49000-0000
300		Black	12× 2 racks, 96 each, blistered	2.304	49000-0006
1100		Black	10× 10 racks, 96 each, cardboard	9.600	49001-0000
1100	+	Black	12× 2 racks, 96 each, blistered	2.304	49001-0106

J.T.BAKER® ROBOTIC TIPS FOR BECKMAN/OLYMPUS

Available in conductive version, with or without filter.

- 140 tips are sealed in one optimized rack
- Ideal for blood analysis
- Working volumes of 20 - 1200 µl


Type `Olympus` tips are compatible to the following platforms:

OLA2500 and AutoMate 1250/2550



J.T.Baker® Conductive tips 1200 µl for Olympus instruments - in Rack

Volume (µ l)	Filter	Tip colour	Packed	Pk	Cat. No.
1200		Black	1 x 8 x 140 each, foiled	1.120	49044-0000
1200	+	Black	1 x 8 x 140 each, foiled	1.120	49044-0100



J.T.Baker® Conductive tips 50 µl for QIAGEN instruments

J.T.BAKER® ROBOTIC TIPS FOR QIAGEN

Available in conductive version, with filter, non-sterile.

- Manufactured under cleanroom conditions (cleanroom class 8)
- State of the art quality control
- Ideal for high-throughput analysis and DNA purification
- Packed in blister

Type `QIAGEN` tips are compatible to the following platforms:

QIAGEN QIAGENity.

Volume (µ l)	Filter	Tip colour	Packed	Pk	Cat. No.
50	+	Black	12x 2 racks, 96 each, blistered	2.304	49018-0106



J.T.Baker® Clear Tips 300 µl for Dynex instruments - blue Sample Tips in Container

J.T.BAKER® ROBOTIC TIPS FOR DYNEX INSTRUMENTS

Available in non-conductive version, without filter, in container rack.

- 108 tips in one optimized container with lid
- Ideal for the clinical diagnostic as well as for drug discovery
- Especially developed for the automated sample handling
- Highest quality specifications for the compatibility and cleanliness
- Blue sample tips for visual identification

Type `Dynex` tips are compatible to the following instruments:

DYNEX DSX and DS2.

Volume (µ l)	Tip colour	Packed	Pk	Cat. No.
300	Blue	1 x 4 racks, 108 each	432	49136-0000
1300	Clear	1 x 4 racks, 108 each	432	49137-0000

Solutions for nucleic acid preparation

DOWNLOAD
YOUR COPY
VWR.COM

Featuring tools needed along the entire workflow:
Sample disruption and homogenisation, nucleic acid
isolation, photometry, centrifugation and storage.



THERMOSHAKE, THERMAL SHAKE TOUCH

The Thermal Shake Touch is designed for applications that require consistent, precise results. With heating and shaking capabilities, this low profile unit uses interchangeable blocks to accommodate a wide variety of tubes and microplates. The LCD touch screen enables faster setting of temperature, speed and time, which can all be viewed at once. Display features on-screen help topics with operational tips. Touch screen is compatible with rubber gloves used in laboratories.

- Easy to use 4,3" colour LCD touch screen allows the user to save, and visibly track progress through the live status bar, for five user-defined programs, each with five individual steps
- Suitable for rapid heating, high speed shaking and a pulse mode feature, ideal for quick vortex applications
- Timer with audible alarm, 1 minute to 99 hours, 59 minutes, heat function will automatically shut off if the unit recognises an internal issue
- USB port can transfer information to a flash drive for data logging, program storage and software updates

Program control capabilities allow user-programmable operation for automated use and memory for five separate, five-step programs. Adjustable temperature ramp rate feature separately defines the heating rates in increments of 0,5 °C/min. Single-point calibration mode for maximum temperature accuracy, the single-point calibration procedure allows the user to calibrate up to six different defined temperatures. Constructed from a high quality heat and chemically resistant polymer, so the housing remains cool to the touch throughout normal operating temperatures. Maximum temperature-limiting function ensures the temperature will not exceed user-defined limits, allowing control of temperature-sensitive samples. A hot top warning illuminates when the temperature reaches 40 °C, and remains on until the unit is cooled below 40 °C.

The unit's enhanced electronics and temperature sensor provide accurate, dependable temperature settings across the operating range of 5 to 35 °C, (maximum 80% relative humidity, non condensing). Applications include cell cultures, DNA, RNA and protein studies.

Supplied with 1,5 ml block (460-0210), rack and cover, additional blocks must be ordered separately. Note that Eppendorf Thermomixer R® blocks are compatible with the VWR® Thermal Shake Touch. Model with NIST traceable certificate is also available, this includes a 3-point NIST traceable calibration. The traceable certificate includes actual calibration measurement data and uncertainty. The calibration laboratory is ISO/IEC 17025 compliant.

*RT = Ambient

** Microtube blocks include a removable rack and cover



Model	Thermal Shake Touch
Heating speed (°C/min)	5 °C/min
Orbit (mm)	3 mm
Speed accuracy (%)	±2%
Speed range (min⁻¹)	300 - 3000 min ⁻¹
Temperature accuracy (°C)	±1 °C (between 20 and 45 °C) ±2 °C (above 45 °C)
Temperature control range (°C)	RT* +4...100 °C
Weight (kg)	4,4 kg
WxDxH (mm)	248x260x132 mm

Description	Pk	Cat. No.
Thermal Shake Touch, EU-plug	1	460-0202
Thermal Shake Touch, UK-plug	1	460-0204
Thermal Shake Touch, CH-plug	1	460-0206
Thermal Shake Touch, NIST certificate, EU-plug	1	460-0203
Thermal Shake Touch, NIST certificate, UK-plug	1	460-0205
Thermal Shake Touch, NIST certificate, CH-plug	1	460-0207

Description	Well size (mm)	For	No. of holes	Depth (mm)	Pk	Cat. No.
Interchangeable blocks for Thermal Shake Touch and Cooling Thermal Shake Touch						
Interchangeable tube block	Ø 7,9	30x0,5 ml tubes**	30	24,6	1	460-0209
Interchangeable tube block	Ø 11,1	24x1,5 ml tubes**	24	35,3	1	460-0210
Interchangeable tube block	Ø 11,5	24x2,0 ml tubes**	24	35,3	1	460-0211
Interchangeable tube block	Ø 12,3	24x12 mm tubes	24	35,3	1	460-0212
Interchangeable tube block	Ø 12,6	24x2,0 ml cryotubes	24	30,5	1	460-0213
Interchangeable tube block	Ø 16,8	9x5 ml Eppendorf tubes	9	49	1	460-0329
Interchangeable tube block	Ø 30,0	4x50 ml conical tubes	4	100,9	1	460-0215
Interchangeable tube block	Ø 17,3	9x15 ml conical tubes	9	104,4	1	460-0214
Interchangeable PCR plate thermal block with lid	Ø 6,4	96x0,2 ml PCR tubes	96	12,7	1	460-0330
Interchangeable 384 well PCR plate thermal block with lid	Ø 4,0	1x384-well PCR microplate	384	8,1	1	460-0331
Interchangeable thermal microplate block with lid	129,5x78,7	1x96-well microplate	1	23	1	460-0208

THERMO SHAKERS, COOLING THERMAL SHAKE TOUCH

The cooling thermal shake touch is designed for applications that require consistent, precise results. With heating, cooling and shaking capabilities, this low profile unit uses interchangeable blocks to accommodate a wide variety of tubes and microplates. The LCD touch screen enables faster setting of temperature, speed and time, which can all be viewed at once. Display features on-screen help topics with operational tips. Touch screen is compatible with rubber gloves used in laboratories.

- Easy to use 4,3" colour LCD touch screen allows the user to save and visibly track progress through the live status bar for five user-defined programs, each with five individual steps
- Suitable for rapid heating, cooling and high speed shaking and a pulse mode feature, ideal for quick vortex applications
- Timer with audible alarm, 1 minute to 99 hours, 59 minutes, heat function will automatically shut off if the unit recognises an internal issue
- USB port can transfer information to a flash drive for data logging, program storage and software updates

Program control capabilities allow user-programmable operation for automated use and memory for five separate, five-step programs. Adjustable temperature ramp rate feature separately defines the heating and cooling rates in increments of 0,5 °C/min. Single-point calibration mode for maximum temperature accuracy, the single-point calibration procedure allows the user to calibrate up to six different defined temperatures. Constructed from a high quality heat and chemical resistant polymer, so the housing remains cool to the touch throughout normal operating temperatures. Maximum temperature-limiting function



ensures the temperature will not exceed user-defined limits, allowing control of temperature-sensitive samples. A hot top warning illuminates when the temperature reaches 40 °C, and remains on until the unit is cooled below 40 °C.

The unit's enhanced electronics and dual-temperature sensors provide accurate, dependable temperature settings across the operating range of 5 to 35 °C, (maximum 80% relative humidity, non condensing). Applications include cell cultures, DNA, RNA, hybridisation and protein studies.

Supplied without blocks, blocks must be ordered separately. Note that Eppendorf Thermomixer R® blocks are compatible with the VWR® Cooling thermal shake touch. Model with NIST traceable certificate is also available, this includes a 3-point NIST traceable calibration. The traceable certificate includes actual calibration measurement data and uncertainty. The calibration laboratory is ISO/IEC 17025 compliant.

** Microtube blocks include a removable rack and cover

Model	Cooling Thermal Shake Touch		
Cooling rate (°C/min)	Above ambient: 2 - 3 Below ambient: 0,5 - 1,0	Above ambient: 2 - 3 °C/min Below ambient: 0,5 - 1,0 °C/min	Above ambient: 2 - 3 Below ambient: 0,5 - 1,0
Heating speed (°C/min)	5		
Orbit (mm)	3		
Speed accuracy (%)	±2		
Speed range (min ⁻¹)	300 - 3000		
Temperature accuracy (°C)	±0,5 (between 20 and 45) ±2 (below 20 and above 45)		
Temperature control range (°C)	From 17 below ambient to 100		
Weight (kg)	4,4		
WxDxH (mm)	248x260x132		

Description	Pk	Cat. No.
Thermo shakers		
Cooling thermal shake touch, EU-plug	1	460-0196
Cooling thermal shake touch, UK-plug	1	460-0198
Cooling thermal shake touch, CH-plug	1	460-0200
Thermo shaker with NIST traceable certificate		
Cooling thermal shake touch, EU-plug, NIST certificate	1	460-0197
Cooling thermal shake touch, UK-plug, NIST certificate	1	460-0199
Cooling thermal shake touch, CH-plug, NIST certificate	1	460-0201

Description	Well size	For	No. of holes	Depth	Pk	Cat. No.
Interchangeable blocks for Thermal Shake Touch and Cooling Thermal Shake Touch						
Interchangeable tube block	Ø 7,9 mm	30x0,5 ml tubes**	30	24,6 mm	1	460-0209
Interchangeable tube block	Ø 11,1 mm	24x1,5 ml tubes**	24	35,3 mm	1	460-0210
Interchangeable tube block	Ø 11,5 mm	24x2,0 ml tubes**	24	35,3 mm	1	460-0211
Interchangeable tube block	Ø 12,3 mm	24x12 mm tubes	24	35,3 mm	1	460-0212
Interchangeable tube block	Ø 12,6 mm	24x2,0 ml cryotubes	24	30,5 mm	1	460-0213
Interchangeable tube block	Ø 16,8 mm	9x5 ml Eppendorf tubes	9	49 mm	1	460-0329
Interchangeable tube block	Ø 30,0 mm	4x50 ml conical tubes	4	100,9 mm	1	460-0215
Interchangeable tube block	Ø 17,3 mm	9x15 ml conical tubes	9	104,4 mm	1	460-0214
Interchangeable PCR plate thermal block with lid	Ø 6,4 mm	96x0,2 ml PCR tubes	96	12,7 mm	1	460-0330
Interchangeable 384 well PCR plate thermal block with lid	Ø 4,0 mm	1x384-well PCR microplate	384	8,1 mm	1	460-0331
Interchangeable thermal microplate block with lid	129,5x78,7 mm	1x96-well microplate	1	23 mm	1	460-0208

Equipment Services

Keep your laboratory in compliance
and in operation



Providing a carefree equipment lifecycle with Avantor Services

- Installation
- Preventive maintenance
- Warranty and repair
- Calibration
- Qualification, validation and certification (IQ/OQ/PQ)
- Equipment management
- Training

[VWR.COM/EIS](https://www.vwr.com/eis)



THERMOSHAKER, THERMAL SHAKE *LITE*

The combination of heating/cooling and shaking makes the Thermal Shake *lite* microtube shaking incubator ideal for many life science research applications in molecular biology, biochemistry and clinical chemistry. Its compact footprint incorporates an intuitive control panel with large multicolour display, allowing users to easily program and view temperature, time and speed settings.

- Choice of eight interchangeable aluminium blocks accommodate PCR plates and tubes ranging from 0,2 to 15 ml
- Fine tune speed control
- Rapid heating and cooling
- Compact footprint

Supplied without blocks, blocks must be ordered separately.

Model	Thermal Shake <i>lite</i>
Heating speed (°C/min)	6,5 °C/min
Orbit (mm)	3 mm
Speed range (min⁻¹)	300 - 1500 min ⁻¹
Temperature accuracy (°C)	±0,5 °C
Temperature control range (°C)	14 °C below ambient to 100 °C
Weight (kg)	8,5 kg
WxDxH (mm)	330x166x240 mm

Description	Pk	Cat. No.
Thermal Shake <i>lite</i>	1	460-0249

Description	For	No. of holes	Pk	Cat. No.
Heating blocks, aluminium, designed for 1,5 ml tubes				
Aluminium block	0,5 ml tubes	54	1	460-0251
Aluminium block	1,5 ml tubes	40	1	460-0250
Aluminium block	0,5 and 1,5 ml tubes	26 + 24	1	460-0255
Aluminium block	2,0 ml tubes	40	1	460-0256
Aluminium block	15 mm tubes	24	1	460-0253
Aluminium block	0,2 ml PCR tubes or plates	96	1	460-0252
Aluminium block	Water bath block (115x73x38 mm)	-	1	460-0254
Aluminium block	96-well ELISA plate	-	1	460-0257

Setting science in motion to create a better world



AUSTRIA

VWR International GmbH
Graumannsgasse 7
1150 Wien
Tel.: +43 1 97 002 0
info.at@vwr.com

BELGIUM

VWR International bv
Researchpark Haasrode 2020
Geldenaaksebaan 464
3001 Leuven
Tel.: +32 (0) 16 385 011
vwr.be@vwr.com

CZECH REPUBLIC

VWR International s. r. o.
Veetee Business Park
Pražská 442
CZ - 281 67 Stříbrná Skalice
Tel.: +420 321 570 321
info.cz@vwr.com

DENMARK

VWR International A/S
Tobaksvejnen 21
2860 Søborg
Tel.: +45 43 86 87 88
info.dk@vwr.com

FINLAND

VWR International Oy
Valimotie 17-19
00380 Helsinki
Tel.: +358 (0) 9 80 45 51
info.fi@vwr.com

FRANCE

VWR International S.A.S.
Le Périgares – Bâtiment B
201, rue Carnot
94126 Fontenay-sous-Bois cedex
Tel.: 0 825 02 30 30* (national)
Tel.: +33 (0) 1 45 14 85 00
(international)
info.fr@vwr.com
* 0,18 TTC/min + prix appel

GERMANY

VWR International GmbH
Hilpertstraße 20a
D - 64295 Darmstadt
Tel.: 0800 702 00 07* (national)
Tel.: +49 (0) 6151 3972 0 (international)
info.de@vwr.com
*Freecall

HUNGARY

VWR International Kft.
Simon László u. 4.
4034 Debrecen
Tel.: +36 52 521130
info.hu@vwr.com

IRELAND

VWR International Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

THE NETHERLANDS

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

NORWAY

VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

PORTUGAL

VWR International - Material de
Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Amadora
Tel.: +351 21 3600 770
info.pt@vwr.com

SPAIN

VWR International EuroLab S.L.U.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
info.es@vwr.com

SWEDEN

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 (0) 8 621 34 00
kundservice.se@vwr.com

SWITZERLAND

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 (0) 44 745 13 13
info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
uksales@vwr.com

CHINA

VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: 400 821 8006
info_china@vwr.com

INDIA

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India - 560058
Tel.: +91 80 28078400
vwr_india@vwr.com

KOREA

VWR International -
17, Daehak 4-ro, Yeongtong-gu,
Suwon-si, Gyeonggi-do
Tel.: +82 31 645 7256
saleskorea@avantorsciences.com

MIDDLE EAST & AFRICA

VWR International FZ-LLC
Office 203, DSP Lab Complex,
Dubai Science Park,
Dubai, United Arab Emirates
Tel.: +971 4 5573271
Info.mea@vwr.com

SINGAPORE

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
sales.sg@vwr.com

GO TO [VWR.COM](https://www.vwr.com) FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS FROM YOUR LOCAL VWR SUPPORT TEAM