












GRAPHITE TUBES FOR AA - SYSTEM

Agilent (Varian)	Description	Original ID	Drawing	Order no. Qty. / VE
	Tube, pyrocoated, with rings recommended standard tube, high quality	63-100012-00		56VA002 10 pcs
	Tube, pyro coated, with rings, Longlife better coating, recommended for over-night-runs	63-100012-EL		56VAEL2 10 pcs.
	Plateau-tube, coated also for use with bone platform	63-100011-00		56VA003 10 pcs.
	Bone-Platform Pyro graphite, inserting in Plateau-tube, not recommended!	63-100013-00		56VA004 10 pcs.
	Tube, uncoated, with rings	63-100015-00		56VA005 10 pcs.
	Tube, coated, with Omega-Platform	63-100026-00		56VA010 10 pcs.
	Plateau-tube, uncoated also for use with bone platform	63-100014-00		56VA033 10 pcs.
	Electrodes / contacts for SpectrAA set of 2 pieces	63-100016-00		56VA006 2 pcs.
	Zeeman-Contacts set of 2 pieces	63-100017-00		56VA007 2 pcs.
	Shroud for SpectrAA 1 piece	63-100018-00		56VA008 1 pc.
	Zeeman Shroud 1 piece	63-100019-00		56VA009 1 pc.
	Zeeman Shroud GTA120	63-100036-00		56VAS1 1 pc.
	Electrodes GTA 120	63-100035-00		56VAS2 2 pcs.

The quality of your AAS graphite tubes is very important to get good results. This is why we only offer tubes made in Germany – OEM quality at the best possible price.

www.hosmotic.com

Hosmotic S.r.l.

Via Raffaele Bosco, 178 – 80069 Vico Equense – NA

P. IVA: 07672391211 – REA: NA-901340

Tel/Fax: 0818798336 - e-mail: info@hosmotic.com